

PDMS Fabrication

Kenneth Schackart¹, Kattika Kaarj¹

¹University of Arizona

[dx.doi.org/10.17504/protocols.io.yadfsa6](https://doi.org/10.17504/protocols.io.yadfsa6)

481b Laboratory

 Kenneth Schackart 

ABSTRACT

This protocol details how to make a small piece of polydimethylsiloxane (PDMS)

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

Lab coat and gloves must be worn at all times

MATERIALS TEXT


Materials

- Delicate Task wipes (Kimwipes®)
- 100 µL pipette and tips
- 10 µL pipette and tips
- 1 mL syringe
- 2 mL tube
- Petri dish lid (35 mm)
- Dessicator in fume hood
- Hot plate
- Aluminum foil
- Blade

Reagents

- Tridecafluoro-1,1,2,2-tetrahydrooctyl trichlorosilane
- PDMS curing agent (silicone elastomer curing agent)
- PDMS polymer (silicone elastomer base)
- 70% ethanol solution

Prepare surface

- 1 Wipe the inside of the petri dish lid with 70% ethanol solution using a Kimwipe® to rid the surface of debris.
- 2 Place petri dish lid(s) and a 10 cm × 10 cm piece of foil in the dessicator, in the fume hood.
- 3 Add  2 µl Tridecafluoro-1,1,2,2-tetrahydrooctyl trichlorosilane to the piece of aluminum foil.

▲ SAFETY INFORMATION


Chemical is highly corrosive, do not inhale or get on skin.

- 4 Dessicate under vacuum for  00:15:00 to vaporize reagent.



Press down on the lid as you apply vacuum. Once vacuum is achieved the lid should be able to be lifted.

Mix polymer

- 5 Add  100 μ l PDMS curing agent and  900 μ l PDMS polymer to 2 mL tube.

- 6 Mix the solution with an inoculation loop until there are small bubbles throughout the solution. This will take approximately  00:05:00

Degas polymer

- 7 Pour the mixed PDMS into petri dish, and try to get all of the polymer out of the tube using a pipette tip or inoculation loop. Evenly spread the polymer in the petri dish by tilting.

- 8 Dessicate under vacuum for  00:30:00 .

Cure polymer

- 9 Place petri dish on hot plate at  120 °C for  00:10:00 covered in aluminum foil.

Harvest polymer

- 10 Carefully peel PDMS from petridish.

- 11 Using a blade, cut out a piece about 1 cm \times 1 cm and set on glass microscope slide. This will be your PDMS surface for water contact angle analysis.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited