

# Immunofluorescence Microscopy Protocol with Methanol Fixed Cells(about S100A6 in SW480 Cells)

Shanshan Feng,Huihong Zhai

## Abstract

**Citation:** Shanshan Feng,Huihong Zhai Immunofluorescence Microscopy Protocol with Methanol Fixed Cells(about S100A6 in SW480 Cells). **protocols.io**

dx.doi.org/10.17504/protocols.io.mqcc5sw

**Published:** 18 Jan 2018

## Protocol

### Sample preparation

#### Step 1.

1. The cells were grown on coverslips for 2 days in DMEM at 37°C and treated with corresponding pretreatment. At the time of fixation, cells should be 50% confluent.
2. Rinse cells briefly in 1X PBS.
3. Fix cells with 4% cold paraformaldehyde for 20 min at room temperature.
4. permeabilized for 10 min with 0.3% Triton X-100 at 37 °C.
5. Rinse three times in 1X PBS, 5 minutes each.
6. Block samples in 5% normal goat serum for 30 min at room temperature.

### Sample staining

#### Step 2.

- 7.Dilute the primary antibody to the recommended concentration/dilution in 5% FBS/PBS overnight at 4 °C.
- 8.Rinse three times in 1X PBS, 5 minutes each.
- 9.Prepare FITC-conjugated secondary antibody in 5% FBS/PBS in the dark for 1 hour at room temperature according to the recommended manufacturer specification data sheet.
- 10.Nuclei were stained with DAPI for 5 min.
- 11.Rinse three times in 1X PBS, 5 minutes each.
- 12.Coverslip with anti-fade mounting medium.
- 13.Seal the edges of the coverslip to the slide to prevent movement of the coverslip while imaging. One sealant that is recommended is a 1:1:1 ratio mixture of vasoline, lanolin and paraffin.
- 14.Immunofluorescence was analyzed with a META-510 Laser Scanning Confocal Imaging System.

**Step 3.**

15. Chamber slides and coverslips.

16. Antibody dilution solution: 5% FBS in 1X PBS