

# Modified total RNA extraction for *Heterosigma akashiwo* using the Qiagen RNeasy kit

Monica Accerbi, Vinay K Nagarajan, Kathryn Coyne, Pamela J Green

## Abstract

This protocol is a modified version of the Qiagen RNeasy assay for isolating total RNA from *Heterosigma akashiwo*. We used this protocol to obtain high quality RNA for downstream analysis such as quantitative PCR.

**Citation:** Monica Accerbi, Vinay K Nagarajan, Kathryn Coyne, Pamela J Green Modified total RNA extraction for *Heterosigma akashiwo* using the Qiagen RNeasy kit. **protocols.io**  
dx.doi.org/10.17504/protocols.io.hipb4dn

**Published:** 30 Sep 2017

## Protocol

### Step 1.

Starting material: 1 ml culture with about 400,000 cells

### Step 2.

To the cells add 1 ml of Buffer RLT with 10  $\mu$ l  $\beta$ -mercaptoethanol

### Step 3.

Centrifuge through a QIAshredder column (750  $\mu$ l at a time) at 11,000 rpm for 1' and collect FT (flow through) in a 5 ml tube

### Step 4.

To the collected FT add 2 ml of 70% (v/v) ethanol and mix at room temperature

### Step 5.

Centrifuge through a RNeasy spin column (750  $\mu$ l at a time) at 11,000 rpm for 1' and discard FT

### Step 6.

Add 750  $\mu$ l Buffer RW1 and centrifuge at 11,000 rpm for 1' and discard FT

### Step 7.

Wash column with 500 µl Buffer RPE, centrifuge at 11,000 rpm for 1' and discard FT

#### **Step 8.**

Repeat step 7

#### **Step 9.**

Centrifuge once more at 11,000 rpm for 1' to dry the membrane and remove residual ethanol

#### **Step 10.**

Elute total RNA with 15 µl DEPC (Diethyl pyrocarbonate) -treated H<sub>2</sub>O

#### **Expected yield**

#### **Step 11.**

Total RNA yield using this method is around 350 ng with 260:280 nm  $\geq$  1.8

#### **For larger volumes of starting material**

#### **Step 12.**

Centrifuge the starting material at 2,000 to 3,000 rpm for 2' at room temperature

Discard supernatant and resuspend pellet in the residual liquid, about 1 ml

Proceed as above, steps 1 to 9. Elute total RNA with 25-30 µl DEPC-treated H<sub>2</sub>O