

Electroporation Protocol (C2986)

New England Biolabs

Abstract


This electroporation protocol is for use with the NEB Turbo Electrocompetent E. coli cells ([C2986](#)). These cells are suitable for high efficiency electroporation and rapid colony growth; they are ideal for DNA library constructions and all cloning purposes.

Citation: New England Biolabs Electroporation Protocol (C2986). [protocols.io](#)

[dx.doi.org/10.17504/protocols.io.crgv3v](https://doi.org/10.17504/protocols.io.crgv3v)

Published: 03 Feb 2015

Materials

 NEB Turbo Electrocompetent E.coli - 6x0.1 ml [C2986K](#) by [New England Biolabs](#)

Protocol

Step 1.

Prepare 17mm x 100mm round-bottom culture tubes (e.g. VWR #60818-667) at room temperature

Step 2.

Place SOC recovery medium in a 37°C water bath

Step 3.

Pre-warm selective plates at 37°C for 1 hour

 [DURATION](#)

01:00:00

Step 4.

Place electroporation cuvettes (1mm) and microcentrifuge tubes on ice

Step 5.

As a positive control for transformation, dilute the control pUC19 by 1:5 to a final concentration of 10 pg/μl using sterile water

[NOTES](#)

New England Biolabs 03 Feb 2015

Heat-denatured ligation reactions can be used for electroporation directly; however, column purification is recommended.

Step 6.

Thaw NEB Turbo Electrocompetent cells on ice (about 10 min) and mix cells by flicking gently

 [DURATION](#)

00:10:00

Step 7.

Transfer **25 μl** of the cells (or the amount specified for the cuvettes) to a chilled microcentrifuge tube

📄 AMOUNT

25 µl Additional info:

Step 8.

Add **1 µl** of the DNA solution

📄 AMOUNT

1 µl Additional info:

Step 9.

Carefully transfer the cell/DNA mix into a chilled cuvette without introducing bubbles and make sure that the cells deposit across the bottom of the cuvette

📄 AMOUNT

1 µl Additional info:

Step 10.

Electroporate using the following conditions for BTX ECM 630 and Bio-Rad GenePulser electroporators: 2.1 kV, 100 Ω, and 25 µF. The typical time constant is 2.6 milliseconds

Step 11.

Immediately add **975 µl** of 37°C SOC to the cuvette

📄 AMOUNT

975 µl Additional info:

Step 12.

Gently mix up and down twice

Step 13.

Transfer to the 17mm x 100mm round-bottom culture tube

Step 14.

Shake vigorously (250 rpm) or rotate at 37°C for 1 hour

🕒 DURATION

01:00:00

Step 15.

Dilute the cells as appropriate then spread 100-200 µl cells onto a pre-warmed selective plate

Step 16.

Incubate plates 8 hours to overnight at 37°C.

🕒 DURATION

08:00:00