

# Immunogold labeling of bacterial cells for transmission electron microscopy (TEM) version 2

Jessica Sacher

## Abstract

**Citation:** Jessica Sacher Immunogold labeling of bacterial cells for transmission electron microscopy (TEM).  
**protocols.io**

dx.doi.org/10.17504/protocols.io.mv5c686

**Published:** 25 May 2018

## Materials

- ✓ 1X PBS (Phosphate-buffered saline ) by Contributed by users
- ✓ 70% Ethanol by Contributed by users
- Paraformaldehyde Aqueous Solution -16% [15700](#) by [Electron Microscopy Sciences](#)
- Anti-Rabbit IgG (whole molecule)-Gold antibody produced in goat [Sigma 0162787](#) by [Sigma](#)
- ✓ 5% BSA in PBST by Contributed by users
- ✓ Formvar-coated 300-mesh electron microscopy grids by Contributed by users

## Protocol

### Step 1.

Cut a strip of Parafilm (20-30 cm long), place flat on bench.

### Step 2.

Prepare one 1-mL tube each of 70% ethanol and mqH2O for washing tweezers between transfers.

### Step 3.

Harvest cells by flooding an overnight plate with 5 mL nutrient broth or buffer, using a spreader to dislodge cells, pipetting into a conical flask and mixing gently by inversion (avoid pipetting up and down, vortexing and centrifuging to preserve flagella).

### Step 4.

Set cells to an OD600 of 1-3 in nutrient broth (1-mL volume), place on ice until use.

### Step 5.

Pipette 100 uL cells onto parafilm, use tweezers to float Formvar-coated copper (300 mesh) grid atop droplet, Formvar-side-down (should do this in duplicate, ie. 2 grids per drop), incubate 45 min at RT (between transferring grids, particularly between transfer of grids containing different samples, swirl tweezers 2-3 seconds in EtOH, then in H2O, then wipe dry with a fresh Kimwipe)

### Step 6.

Wash grids 1x by transferring to a 100-uL drop of PBS, incubating 3-5 min.

### Step 7.

Incubate 20 min in 2.5% paraformaldehyde/PBS to kill cells and fix them to grids.

**Step 8.**

Wash 3x in PBS (100 uL drops, 3-5 min each).

**Step 9.**

Block in 5% BSA/PBST (100 uL, 35 min).

**Step 10.**

\*If doing Far Western blot, ie. testing protein binding to cells, incubate in desired protein diluted 1/25 in blocking solution, 45 min, 100 uL drop

**Step 11.**

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

**Step 12.**

Incubate in primary antibody (anti-Gp047, rabbit) at 1/50 (diluted in blocking solution), 45 min, 100 uL.

**Step 13.**

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

**Step 14.**

Incubate in secondary antibody (goat anti-rabbit IgG-gold) at 1/50 (diluted in blocking solution), 45 min, 100 uL.

**Step 15.**

Wash 3x in blocking solution (100 uL drops, 3-5 min each).

**Step 16.**

Wash 3x in PBS (100 uL drops, 3-5 min each).

**Step 17.**

Wash 3x in water (100 uL drops, 3-5 min each) to remove salts, which can form crystals under EM.

**Step 18.**

Dry completely on Whatman filter paper in a Petri dish (at least 20 min, preferably overnight). OK to store grids in this way at RT for several days prior to TEM, but ideally visualize within 1 day.

**Step 19.**

Analyze grids via TEM.