



Oct 06, 2019

Preparation of primary chicken embryo liver (CEL) cells

Norfitriah Mohamed Sohaimi¹, Mohd Hair Bejo¹, Abdul Rahman Omar¹, Aini Ideris¹, Nurulfiza Mat Isa²

¹Faculty of Veterinary Medicine, Universiti Putra Malaysia, Serdang, Selangor, Malaysia., ²Faculty of Biotehnology and Biomolecular Sciences, Universiti Putra Malaysia, Serdang, Selangor, Malaysia.

dx.doi.org/10.17504/protocols.io.7zchp2w 1 Works for me

Norfitriah Mohamed Sohaimi 🕢

MATERIALS

NAME ~	CATALOG #	VENDOR V
Dulbecco's Modified Eagle's Medium (DMEM)	D5796	Sigma Aldrich
Fetal Bovine Serum	10270106	Gibco - Thermo Fischer
Gibco Penicillin-Streptomycin (10,000 U/mL) (Pen/Strep)	15-140-122	Fisher Scientific
0.1M Phosphate Buffered Saline pH 7.4		
Trypsin-EDTA (0.25%), phenol red	25200056	Thermo Fisher

MATERIALS TEXT

Primary CEL cell was obtained from liver embryos of 13 to 15 days old SPF embryonated chicken eggs. Liver was harvested aseptically using sterile forceps and washed twice with sterile phosphate buffered saline (PBS, pH 7.4, 0.1M). The liver tissue was minced and trypsinized gently with 0.25% Trypsin-EDTA solution for 10 minutes. The suspension was passed through muslin cloth and centrifuged at 96 x g for 5 minutes to obtain cell pellet. Trypsin was discarded and cell pellet was resuspended with fresh Dulbecco's Modification Eagle Medium (DMEM), enriched with 10% fetal bovine serum (FBS) and 1% penicillin-streptomycin antibiotic. Cell concentration was counted and adjusted to 5×10^6 cells/mL. The cell suspension (5mL) was seeded into new 25cm^2 cell culture flasks and was kept under controlled atmosphere at 5% CO₂ incubator with 85%-90% humidity until confluent monolayer formed [Soumyalekshmi et al., 2014].

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited