



May 14,
2019

Working

Determination of florfenicol and florfenicol amine in fish plasma (*Salmo salar*) through HPLC MS/MS [↗](#)

PLOS One

Betty San Martín¹, Lisette Lapierre¹, Marcela Fresno¹, Javiera Cornejo¹

¹Universidad de Chile

[dx.doi.org/10.17504/protocols.io.zhdf326](https://doi.org/10.17504/protocols.io.zhdf326)

Lisette Lapierre

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215174>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Martín BS, Fresno M, Cornejo J, Godoy M, Ibarra R, Vidal R, Araneda M, Anadón A, Lapierre L (2019) Optimization of florfenicol dose against *Piscirickettsia salmonis* in *Salmo salar* through PK/PD studies. PLoS ONE 14(5): e0215174. doi: [10.1371/journal.pone.0215174](https://doi.org/10.1371/journal.pone.0215174)

STEPS MATERIALS

NAME	CATALOG #	VENDOR	CAS NUMBER	RRID
Water				
Acetone				
Dichloromethane	DC3600.SIZE.1L	Bio Basic Inc.		
Water				
Acetoacetic acid				
Water				
Acetoacetic acid				
Water				
Methanol				

Sample preparation

1 Extraction was carried out according to the method described by Li *et al.* (2006), with modifications.

2 250 µl of each plasma sample were placed in 15 ml Falcon tubes and 500 µl



Water

and 2 ml



Acetone

were added.

- 3 The tubes were shaken for ⌚ 00:10:00 , sonicated for ⌚ 00:10:00 and centrifuged at 5,000 rpm for ⌚ 00:05:00

- 4 Subsequently, the resulting supernatant was transferred to another 🧴 15 ml Falcon tube containing 🧴 2 ml



Dichloromethane

by Bio Basic Inc.

Catalog #: DC3600.SIZE.1L

, and then shaken, sonicated and centrifuged at 5,000 rpm for the same amount of time as before

- 5 Then, the upper phase (



Water

) was discarded and the remaining phase evaporated at 🔥 40 °C to 🔥 50 °C under nitrogen flow.

- 6 Once completely dry, 🧴 1 ml of mobile phase and 🧴 1 ml hexane were used to reconstitute each sample, which then was shaken, sonicated and centrifuged at 5,000 rpm.

- 7 Finally, 🧴 300 µl of the lower phase were transferred to a vial to be injected into the

- 8 The limit of detection (LOD) and the limit of quantification (LOQ) established for both analytes in plasma were 🧴 0.1 Parts per Million (PPM) and 🧴 0.2 Parts per Million (PPM) , respectively.


Instrumental Analysis

- 9 For instrumental analysis, an LC (Agilent, 1290 infinity series) coupled to a triple quadrupole mass spectrometer (API 5500, ABSCIEX) was used.

10 A Synergi™ 4 µm fusion RP 80Å 50 x 2.0 mm analytic column was used.

11 Analyst 1.6.3 and Multiquant 3.0 software was used for equipment management and integration, respectively.


12 Chromatographic separation was performed through a mobile phase using solvent A: 0.1% of

 Acetoacetic acid

in

 Water

; and a mobile phase solvent B: 0.1% of



 Acetoacetic acid

in

 Water

/

 Methanol

1:9 ratio; with a gradient flow of  350 µl min⁻¹ and a gradient elution from 32% up to 68% solvent A in  00:03:00 of 35% phase solvent A, and 75% phase solvent B.

13 Injection volume was  2 µl; column oven temperature was set at  37 °C.

14 Parameters of the MS/MS detector

Ionization	Electrospray (ESI)
Scan type	MRM
Source temperature (TEM)	550°C
Nebulizer (GS1)	60 psi
Turbo ion (GS2)	80 psi
Curtain gas (CUR)	20 psi
Collision gas (CAD)	10 psi
Ion-spray voltage (IS)	4500 V

15 Monitored ion masses

Analyte	Precursor ion (Q1 mass) (Da)	Fragment ion (Q3 mass) (Da)	Time (ms)	EP (V)	CE (V)	CXP (V)
FF1	356.0	336.0	100.0	-5,000	-15,000	-8,000
FF2	356.0	185.0	100.0	-5,000	-17,000	-12,000
FFA1	248.0	230.0	200.0	5,000	22,000	25,000
FFA2	248.0	130.0	200.0	2,000	30,000	10,000
CAF-d5 (IS)	326.0	157.0	100.0	10,000	-25,000	-20,000

Validation of analytical method

- 16 Prior to determining concentrations in plasma, the analytical method was validated by HPLC MS/MS, according to instructions from the European Commission Decision 2002/657/EC (2002).
- 17 Precursor ions and the two product ions were identified for florfenicol and florfenicol amine, respectively.
- 18 Values for essential parameters were estimated for the validation of the analytical method on plasma: specificity, recovery, repeatability, intralaboratory reproducibility, decision limit (CC α), detection capability (CC β), and linearity.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited