



Mar 23,  
2019

In devel.

## 1. RNA isolation for tissue

[Chin Yee Tan](#)<sup>1</sup>

<sup>1</sup>Duke University

[dx.doi.org/10.17504/protocols.io.zfyf3pw](https://doi.org/10.17504/protocols.io.zfyf3pw)

 [Chin Yee Tan](#) 

### PROTOCOL STATUS

#### In development

We are still developing and optimizing this protocol

- 1 Homogenization in TriZol  
Add 1mL of Trizol reagent per 30mg of tissue and homogenize using handheld homogenizer
- 2 Incubate at RT for 5 mins to allow nucleoprotein complexes to dissociate
- 3 Add 1/5 the volume of trizol of chloroform carefully, and vortex to mix well
- 4 Spin down at max speed, chilled centrifuge for 15 minutes
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube  
The interphase and bottom organic phase can be saved for DNA and protein respectively
- 6 To the aqueous phase, add 500 microliters of 100% isopropanol, mix by inversion and incubate at -20C for a minimum for 2hrs
- 7 Spin down at max speed for 30 minutes to precipitate RNA
- 8 Remove supernatant, and add 1mL of 75% Ethanol to wash the pellet
- 9 Spin down at max speed for 15 minutes and remove supernatant
- 10 Resuspend pellet in appropriate volume of nuclease free water



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited