# **Pseudoalteromonas Media Recipes**

#### **Matthew Sullivan**

#### **Abstract**

Citation: Matthew Sullivan Pseudoalteromonas Media Recipes. protocols.io

dx.doi.org/10.17504/protocols.io.de23gd

Published: 02 Feb 2016

#### **Guidelines**

### Zobell plates, 100% nutrient (1 L)

26 g sea salts 1 g yeast extract 5 g proteose peptone no. 2 fill to 1 L  $qH_2O$ 12 g agar

Grow P.a. in PZM (=50% nutrient)

#### Zobell plates, 20% nutrient (1 L)

26 g sea salts
0.2 g yeast extract
1 g proteose peptone no. 2
fill to 1 L qH<sub>2</sub>O
12 g agar

Autoclave 1 L for 50 min., pour 10 ml/plate

### Zobell liquid, (1 L)

26 g sea salts 1 g yeast extract 5 g proteose peptone no. 2 fill to 1 L qH<sub>2</sub>O

#### Top Zobell agar (; 1 L)

26 g sea salts 1 g yeast extract 5 g proteose peptone no. 2 fill to 1 L qH₂O 6 g agar

#### MSM (1 L)

100 mM NaCl (5.8 g) 81.2 mM MgSO<sub>4</sub>·7H<sub>2</sub>O (20 g) 50 mM Tris base (3.02 g) OR 50 mL 1 M Tris (pH 7.6) fill to 1 L  $qH_2O$  adjust to pH 7-8 if necessary

Artificial seawater, 26 psu (ASW; 1 L) 26 g sea salts fill to 1 L  $qH_2O$ 

## **Protocol**