

# Transformation of Bacterial Cultures Using Hexamine Cobalt Chloride

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova Transformation of Bacterial Cultures Using Hexamine Cobalt Chloride. **protocols.io**

dx.doi.org/10.17504/protocols.io.eq7bdzn

**Published:** 13 Jun 2016

## Guidelines

### Materials:

#### 1) SOB medium

- 1.0% Bacto-tryptone
- 0.5% Bacto-yeast extract
- 10.0 mM NaCl
- 2.5 mM KCl
- 10.0 mM MgCl<sub>2</sub>
- 10.0 mM MgSO<sub>4</sub>
- Prepare the MgCl<sub>2</sub> and MgSO<sub>4</sub> as 1 M stock solutions and autoclave separately.
- Add the MgCl<sub>2</sub> and MgSO<sub>4</sub> after sterilization of the remainder of the components.

#### 2) SOC medium

- SOB media, supplemented with 20mM glucose

#### 3) TFB buffer

- 10 mM K-MES, pH 6.2, 100 mM KCL, 45 mM MnCl<sub>2</sub>·4H<sub>2</sub>O, 10 mM CaCl<sub>2</sub>·2H<sub>2</sub>O, 3 mM NaCoCl<sub>3</sub>
- Weigh out the components. Add the MES to d-H<sub>2</sub>O and adjust the pH with KOH.
- Add the remaining components (in order), waiting until one component is in solution before adding the next. Adjust the volume to the final volume. Filter sterilize. Store frozen at -20°C in 15 mL aliquots.

#### 4) DTT solution

- 2.25 M DTT, 40 mM KOAc, pH 6.0. Filter sterilize, Store frozen at -20°C.

#### 5) DMSO (Dimethylsulfoxide)

## Reference

D. Hanahan. (1983). Studies on Transformation of *Escherichia coli* with plasmids. Journal of Molecular Biology **166**: 557-580.

## Protocol

### Step 1.

Prepare the SOB medium.

✓ [PROTOCOL](#)

#### . [SOB medium](#)

CONTACT: [Irina Agarkova](#)

#### Step 1.1.

Prepare the  $\text{MgCl}_2$  and  $\text{MgSO}_4$  as 1 M stock solutions and autoclave separately.

#### Step 1.2.

Add the  $\text{MgCl}_2$  and  $\text{MgSO}_4$  after sterilization of the remainder of the components.

### Step 2.

Prepare the SOC medium by using the SOB media supplemented with 20mM glucose.

### Step 3.

Prepare the TFB buffer.

✓ [PROTOCOL](#)

#### . [TFB buffer](#)

CONTACT: [Irina Agarkova](#)

#### Step 3.1.

Weigh out the components.

#### Step 3.2.

Add the MES to d- $\text{H}_2\text{O}$  and adjust the pH with KOH.

#### Step 3.3.

Add the remaining components (in order), waiting until one component is in solution before adding the next.

#### Step 3.4.

Adjust the volume to the final volume.

#### Step 3.5.

Filter sterilize.

#### Step 3.6.

Store frozen at  $-20^\circ\text{C}$  in 15 ml aliquots.

### Step 4.

Prepare the DTT solution.

## **PROTOCOL**

### **. DTT solution**

CONTACT: [Irina Agarkova](#)

#### **Step 4.1.**

Add components together and adjust pH as necessary.

#### **Step 4.2.**

Filter sterilize.

#### **Step 4.3.**

Store frozen at -20°C.

#### **Step 5.**

Grow 5 mL of the host cells overnight in SOB media at 37°C.

 **DURATION**

18:00:00

#### **Step 6.**

Inoculate 40 mL of SOB medium with 0.8 mL of the overnight culture.

#### **Step 7.**

Grow to an  $A_{550}$  of 0.45-0.55 at 37°C (approximately 3-4 hours).

 **DURATION**

04:00:00

#### **Step 8.**

Centrifuge the cells in the Sorvall SS34 rotor at 5,000 rpm, 5 min, 4°C.

 **DURATION**

00:05:00

#### **Step 9.**

Discard the supernatant.

#### **Step 10.**

Resuspend the pellet with 12.5 mL of the TFB solution.

#### **Step 11.**

Hold the remaining 2.5 mL of TFB for use later.

#### **Step 12.**

Chill the cells on ice for 15 min.

 **DURATION**

00:15:00

#### **Step 13.**

Centrifuge the cells in the Sorvall SS34 rotor at 5,000 rpm, 5 min, 4°C.

 **DURATION**

00:05:00

#### **Step 14.**

Discard the supernatant.

#### **Step 15.**

Resuspend the pellet with 2.4 mL of TFB solution.

**Step 16.**

Add DMSO to 3.5% (84µL), mix and chill on ice for 5 min.

 DURATION

00:05:00

**Step 17.**

Add DTT solution to 75 mM (84 µL), mix and chill on ice for 10 min.

 DURATION

00:10:00

**Step 18.**

Add an equal volume of DMSO as before (84 µL), mix and chill on ice for 5 min. The cells are now "competent".

 DURATION

00:05:00

**Step 19.**

Pipet 21 µL competent cells per prechilled microfuge tube.

 NOTES

**Irina Agarkova** 14 Apr 2016

One tube will be spread on one plate.

**Step 20.**

Add the DNA (in as small a volume as possible, 1-2 µL/tube), mix and chill on ice for 30 min.

 DURATION

00:30:00

**Step 21.**

Heat pulse the tubes at 42°C for 3 min.

 DURATION

00:05:00

**Step 22.**

Then chill on ice for 2 min.

 DURATION

00:02:00

**Step 23.**

Add 80 µL of SOC medium per tube and incubate the tubes at 37°C for 60 min.

 DURATION

01:00:00

**Step 24.**

Spread 100 µL onto each plate.

**Step 25.**

Incubate the plates at 37°C overnight.

 DURATION

18:00:00