



Detection and quantification of *Candida* spp. from subgingival areas [↗](#)

PLOS One

[Sanja Petrovic](#)¹, Milena Radunovic¹, Ana Pucar¹

¹University of Belgrade

[dx.doi.org/10.17504/protocols.io.wdefa3e](https://doi.org/10.17504/protocols.io.wdefa3e)

[Sanja Petrovic](#)

ABSTRACT

Detection and quantification of *Candida* spp. from subgingival areas using two methods- sterile periodontal curette and sterile paper points.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0210527>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Petrovic SM, Radunovic M, Barac M, Pficer JK, Pavlica D, Arsenijevic VA, Pucar A (2019) Subgingival areas as potential reservoirs of different *Candida* spp in type 2 diabetes patients and healthy subjects. PLoS ONE 14(1): e0210527. doi:

[10.1371/journal.pone.0210527](https://doi.org/10.1371/journal.pone.0210527)

PROTOCOL STATUS

Working

MATERIALS

NAME	CATALOG #	VENDOR
Disposable gloves, nitrile		
Periodontal probe	PCPNT2W	
G. Hartzell & Son Model # S4L/4RHSS; Heavy Duty Columbia Curette One End # S4L One End # 4RHSS	126-S4L/4RHSS	
Tab Top Microcentrifuge Tube Sorenson SioScience Inc. 1.5ml	16479	
Sabouraud Dextrose Broth	M033	
Sabouraud dextrose agar	MH063-100G	
Candida differential agar	M1297A-100G	

SAFETY WARNINGS

1. The first day- clinical periodontal examination (full mouth examination at six sites per tooth) in order to diagnose chronic periodontitis and detect the deepest periodontal pocket (with the highest PPD value).
2. -Record plaque index (Silness-Löe), bleeding on probing (BOP), periodontal pocket depth (PPD) and clinical attachment loss (CAL). Samples should not be obtained the first day, to avoid contamination by blood.
3. 2. The second day- sampling. Isolate sampling tooth with cotton rolls.

- 4 -Remove the supragingival plaque by a curette and sterile gauze.
- 5 -Place two sterile paper points (#30) into the pocket/sulcus until a mild resistance appears and keep for 30 seconds. Paper points contaminated by blood should be excluded.
- 6 -Collect paper points and place into 1mL of Sabouraud dextrose broth.
- 7 - Use sterile curette to obtain samples of the complete subgingival biofilm and inoculated in separate sterile plastic tubes containing 1 mL of Sabouraud dextrose broth.-
- 8 - Send samples to microbiologic laboratory.
- 9 3. The second day -microbiological protocols. Vortex samples for 60 seconds.
- 10 - Streak 20µl of suspended broth on SDA and ChromAgar using sterile plastic micro pipette.
- 11 -Incubate at 25°C for 48 hours
- 12 4. The fourth day. Count the number of Colony Forming Units per sample.
- 13 -Multiply the obtained number by 50 in order to obtain CFU/ml.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited