

# Northern blotting with DIG Northern Starter Kit (Roche®)

Lutz Berwanger

## Abstract

This script is used for northern blotting with the Roche DIG Northern Starter Kit.

**Citation:** Lutz Berwanger Northern blotting with DIG Northern Starter Kit (Roche®). **protocols.io**

dx.doi.org/10.17504/protocols.io.guhbw6

**Published:** 27 Mar 2017

## Protocol

### Preparation of needed buffers

#### Step 1.

It is recommended to prepare all needed buffers before start northern blotting.

List of Buffers:

Solution	Preparation	Storage	Use
DIG Easy Hyb Granules	Reconstitue granules adding 64 mL DMPC treated water, dissolve by stirring at 37° C	at +15° C to 25° C for 1 month	As prehybridization and hybridization solution
Washing buffer	0.1 M Maleic acid, 0.15 M NaCl; adjust with NaOH (solid) to pH 7.5 (20° C); 0.3% (v/v) Tween 20	at +15° C to 25° C stable	Removal of unspecific bound antibody
Maleic acid buffer	0.1 M Maleic acid, 0.15 M NaCl; adjust with NaOH (solid) to pH 7.5 (20° C)	at +15° C to 25° C stable	Dilution of Blocking solution
Detection buffer	0.1 Tris-HCl, 0.1 M NaCl, adjustment of pH to 9.5 (20° C)	at +15° C to 25° C stable	Adjustment of pH to 9.5
Blocking solution	Prepare a 10x stock solution	at 8° C stable	Blocking of unspecific binding sites on the membrane

### Preparation of hybridization probe

#### Step 2.

Producing **5S RNA** hybridization probe:

Primer:

TCTTGGTGTCTTTAGCGTCATGGAAC (#460\_5S\_6803\_fw)

TAATACGACTCACTATAGGGACTTGGCATCGGACATTGTGC (#461\_T7\_5S\_6803\_rev)

PCR reaction mix:

5x Q5 Reaction buffer	10 µL
10 mM dNTPs	1 µL
10 µM forward Primer	2.5 µL
10 µM revers Primer	2.5 µL
template DNA	1 µL
Q5 polymerase	0.5 µL
H <sub>2</sub> O	ad 50 µL

PCR setup:	Anneal at 71°C		
	1	98°C	30 sec
Loop 30x			
	2	98°C	5 sec
	3	71°C	15 sec
	4	72°C	15 sec
end of Loop			
	5	72°C	2 min
	6	4°C	10 sec

## ⓘ NOTES

**Lutz Berwanger** 20 Dec 2016

5S RNA hybridization probe is an example. PCR recepies vary of course.

Preparation of denaturing 8% PAA gels

### Step 3.

Clean glassplates, spacer and combs with 70% EtOH and RNase Away. Assemble gelchamber in casting frame.

Gel mixture for 10 mL:

8% PAA gel:

40% PAA (19:1) Rotipage (Roth)	2 mL
Urea	5 g
10x TBE	1 mL
H <sub>2</sub> O	ad 10 mL
10x APS	100 µL
TEMED	10 µL

Pour gel and let polymerize at least 2 h at roomtemperature.

 **DURATION**  
02:00:00

In-vitro-transcription of hybridization probe

#### Step 4.

In-vitro-transcription mixture:

PCR product (0.1-1 µg)	14 µL
T7-Flash 10x reaction buffer	2 µL
DIG RNA labeling mix	2 µL
100 nM DTT	2 µL
RiboGuard RNase Inhibitor	0.5 µL
T7-Flash Enzyme Solution	2 µL

Incubate at 37°C at least 2 h.

 **DURATION**  
02:00:00

 **NOTES**

**Lutz Berwanger** 20 Dec 2016

Start in-vitro-transcription when assembling gel and stop it, when start hybridization step.

RNA preparation and gel loading

#### Step 5.

Prepare RNA-probes:

- Use 3 µg of previously extracted RNA (see mentioned protocol) and aquibrilate probes on same level with DEPC-treated H<sub>2</sub>O and add 2x RNA loadingbuffer (Thermofisher).
- Prepare mixture of 4 µL high range ladder and 2x RNA loadingbuffer (Thermofisher).
- Denature probes and ladder for **5 min at 95°C** and rapidly cool down on ice.
- Load gel and let it run at **25 mA in 1% TBE buffer**.

## **PROTOCOL**

### **RNA extraction with PGTX**

CONTACT: [Lutz Berwanger](#)

## **NOTES**

**Lutz Berwanger** 20 Dec 2016

It is very important to wash pockets of the gel before loading the gel!

### **Step 5.1.**

Centrifuge for 3 minutes at 4°C and 4.000 g.

## **DURATION**

00:03:00

### **Step 5.2.**

Resuspend the cell pellet in 1 ml of PGTX. Freeze in liquid nitrogen and store at -20 °C.

## **SAFETY INFORMATION**

**Wear goggles, a lab coat and gloves when dealing with PGTX and liquid nitrogen.**

### **Step 5.3.**

Incubate for 5 min at 95 °C, shaking at 250 rpm in Thermomixer (Eppendorf)

## **DURATION**

00:05:00

### **Step 5.4.**

Rapidly chill 5 min on ice.

## **DURATION**

00:05:00

### **Step 5.5.**

Add 700 µl Chloroform.

## **REAGENTS**

✓ Chloroform by Contributed by users

#### SAFETY INFORMATION

**Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.**

#### **Step 5.6.**

Let the samples incubate for 10 min at room temperature in a Thermomixer. Vortex from time to time.

#### DURATION

00:10:00

#### **Step 5.7.**

Centrifuge for 15 min at 14.000 g, 4 °C. Transfer the upper aqueous phase (**contains RNA**) to a fresh reaction tube and add the same volume (450 µL) of Aqua- P/C/I (Phenol/Chloroform/Isoamylalcohol).

#### DURATION

00:15:00

#### SAFETY INFORMATION

**Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.**

#### **Step 5.8.**

Thoroughly mix by vortexing. Centrifuge for 15 min at high speed. Transfer the upper aqueous phase to a 1.5 mL reaction tube.

#### DURATION

00:15:00

#### SAFETY INFORMATION

**Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.**

#### **Step 5.9.**

Add 1 volume of isopropanol.

#### REAGENTS

✓ Isopropanol by Contributed by users

#### **Step 5.10.**

Mix and incubate for at least 30 min at -20 °C. (Can be left overnight.)

#### DURATION

00:30:00

#### **Step 5.11.**

Centrifuge at 14.000 g for 30 min.

#### DURATION

00:30:00

#### **Step 5.12.**

Wash with 75% chilled ethanol. Avoid resuspending the pellet.



#### REAGENTS

Ethyl alcohol, Pure 200 proof, for molecular biology [E7023](#) by [Sigma Aldrich](#)

#### Step 5.13.

Discard supernatant.

#### Step 5.14.

Centrifuge at 14.000 g.



#### DURATION

00:05:00

#### Step 5.15.

Repeat washing step with 75% chilled ethanol.



#### NOTES

**Lutz Berwanger** 19 Dec 2016

Remove excess of EtOH by using a pipet.

#### Step 5.16.

Air dry pellet at RT. Do not overdry!



#### DURATION

00:10:00

#### Step 5.17.

Resuspend the pellet with 30 µL volume of ddH<sub>2</sub>O.



#### NOTES

**Lutz Berwanger** 19 Dec 2016

Usually 40 µL DECP-treated H<sub>2</sub>O

#### Step 5.18.

Fill 50 mL Falcon tube with ice. Fill with bacterial liquid culture, up to a volume of 45 mL.

#### Step 5.19.

Discard supernatant. Resuspend cell pellet in residual water (1 mL). Transfer to 2 mL 'safe lock' tube.

#### Step 5.20.

## DURATION

00:00:15

### Semi-Dry blotting with Trans-Blot SD

#### Step 6.

Preparation of Semidry blotting sandwich:

- Prepare the transfer buffer (1x TBE Buffer)
- Cut the membrane to the dimensions of the gel. Wet the membrane in transfer buffer.
- Cut filter paper ( 8 peaces of thick filter paper) to the dimensions of the gel. Wet filter paper in transfer buffer.
- Place 4 of the pre-soaked sheets of thick filter paper onto the platinum anode. Roll a pipet or test tube over the surface of the filter paper to exclude all air bubbles.
- Place the pre-wetted blotting media on top of the filter paper.
- Carefully place the gel on top of the transfer membrane, aligning the gel on the center of the membrane.
- Place the other pre-soaked sheets of filter paper onto the gel. Exclude air bubbles as before.
- Carefully place the cathode onto the stack.

Blot with 2 mA/ cm<sup>2</sup> membrane for 45 min.

## DURATION

00:45:00

### UV-crosslinking

#### Step 7.

UV-crosslink the wet membrane without prior washing. After crosslinking, membrane can be stored dry at +8°C to +25°C. If you want to go ahead use the membrane immediatley for prehybridization.

### Prehybridization

#### Step 8.

- Prewarm DIG Easy Hyb solution to hybridization temperature (62° C).
- Prehybridize membrane with RNA side facing inwards with DIG Easy Hyb for 30 min with gentle agigation in hybridization flask.

## DURATION

00:30:00

## NOTES

**Lutz Berwanger** 06 Jan 2017

Hybridization temperature differs from one probe to each other.

### Hybridization

## Step 9.

- Denature DIG-labeled RNA probe from **Step 4** by boiling for **5 min** and rapidly cooling in ice/water.
- Add denatured DIG-labeled RNA probe to hybridization flask and incubate for **6 h** or **over night** with gentle agitation.

Hybridization is followed by stringency washes:

- Wash **2 x 5 min** in **2x SSC, 0.1% SDS** at room temperature under constant agitation.
- Wash **2 x 15 min** in **0.1x SSC, 0.1% SDS** at **62° C** under constant agitation.

Immunological detect

## Step 10.

Stripping and reprobing of RNA blots

## Step 11.