

# Sueoka's High Salt Medium

Steven Burgess

## Abstract

Sueoka's High Salt Medium is used for photoautotrophic growth of the green alga *Chlamydomonas reinhardtii*. The medium is also known as HS or HSM.

Source:

Sueoka, N. (1960) *Proc. Natl. Acad. Sci. USA* **46**, 83-91.

<http://www.chlamycollection.org/Sueoka.html>

**Citation:** Steven Burgess Sueoka's High Salt Medium. **protocols.io**

[dx.doi.org/10.17504/protocols.io.fdebi3e](https://doi.org/10.17504/protocols.io.fdebi3e)

**Published:** 23 Jul 2016

## Protocol

### Step 1.

Salt Solution

 **AMOUNT**

5 ml Additional info:

 **PROTOCOL**

. [Beijerinck's Solution](#)

CONTACT: [Steven Burgess](#)

#### Step 1.1.

NH<sub>4</sub>Cl

 **AMOUNT**

100 g Additional info:

#### Step 1.2.

CaCl<sub>2</sub> · 2H<sub>2</sub>O

 **AMOUNT**

2 g Additional info:

### Step 1.3.

$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

 [AMOUNT](#)

4 g Additional info:

### Step 1.4.

Distilled  $\text{H}_2\text{O}$

## Step 2.

Phosphate solution

 [AMOUNT](#)

5 ml Additional info:

 [PROTOCOL](#)

. [Phosphate solution](#)

CONTACT: [Steven Burgess](#)

### Step 2.1.

Potassium hydrogen phosphate

$\text{K}_2\text{HPO}_4$

 [AMOUNT](#)

29 g Additional info:

### Step 2.2.

Potassium dihydrogen phosphate

$\text{KH}_2\text{PO}_4$

 [AMOUNT](#)

14 g Additional info:

## Step 3.

Hutner's Trace Elements

 [AMOUNT](#)

1 ml Additional info:

## Step 4.

Distilled  $\text{H}_2\text{O}$

 [NOTES](#)

**Steven Burgess** 21 Jul 2016

For solid medium, add 15 g agar per liter.