

Staining Labyrinthulomycetes with Propidium Iodide

This protocol was modified by Sabrina Geraci-Yee and is from: Pozarowski P, Darzynkiewicz Z. 2004. Analysis of cell cycle by flow cytometry. Methods Mol Bio. 281:301-311.

Abstract

Sabrina Geraci-Yee modified this protocol from: Pozarowski P., Darzynkiewicz Z. 2004. Analysis of cell cycle by flow cytometry. Methods Mol. Bio. 281:301-311.

Citation: This protocol was modified by Sabrina Geraci-Yee and is from: Pozarowski P, Darzynkiewicz Z. 2004. Analysis of cell cycle by flow cytometry. Methods Mol Bio. 281:301-311. Staining Labyrinthulomycetes with Propidium Iodide. protocols.io

dx.doi.org/10.17504/protocols.io.hfjb3kn

Published: 30 Mar 2017

Protocol

Step 1.

📌 NOTES

Laura Halligan 30 Mar 2017

ATCC Medium: 790 By+ Medium

Yeast Extract.....1.0 g
Peptone.....1.0 g
D+-Glucose.....5.0 g
Seawater.....1000 ml

<https://www.atcc.org/~media/920FDAC93FF84B79851C29FBB8049862.ashx>

'Artemia Soup' medium

1.25 g ground, freeze-dried brine shrimp
autoclaved in 1 liter artificial seawater

Step 2.

Re-suspend cell pellet in 1.5 ml 70% EtOH to fix the cells. Store at 4°C for at least 1 hour (and up to several weeks).

Step 3.

Prepare PI staining solution (Keep propidium iodide stain in dark) :

- 0.1% (v/v) Triton X-100 → 1 ul per 1 ml stain
- 10 ug/mL Propidium iodide (PI) → 1 ul per 1 ml stain
- 100 ug/mL DNase-free RNase A → 10 ul per 1 ml stain
- Bring up to 1 ml by adding 1X PBS → 988 ul per 1 ml stain

Step 4.

Centrifuge ethanol-suspended cells at 10,000 x g for 5 minutes, remove ethanol thoroughly. Pour off and blot the top of the tube. If necessary, pipette out remaining ethanol without disturbing pellet.

Step 5.

Re-suspend the cell pellet in 1 ml of PI staining solution; keep in dark at room temp for 30 min or at 37°C for 10 min.

Step 6.

View cells under fluorescence (you may also want to concentrate the cells on a filter). Excitation maximum is 535 nm and emission maximum is 617 nm.