

# Euplotes crassus transformation using FuGene HD Transfection Reagent as vehicle Version 2

RACHELE CESARONI

## Abstract

**Citation:** RACHELE CESARONI Euplotes crassus transformation using FuGene HD Transfection Reagent as vehicle. **protocols.io**

dx.doi.org/10.17504/protocols.io.g3wbype

**Published:** 30 Jan 2017

## Protocol

### Step 1.

$2 \times 10^4$  well-fed cells (we used E.coli as the only food source) were collected by centrifugation (400 rcf for 3 minutes).

### Step 2.

Cells were washed twice with artificial sea water (see attachment for the recipe) and once with 500 mM sorbitol, 0.5 mM Tris-HCl pH 7.0 (400 rcf for 3 minutes each time), in which they were also resuspended in 50  $\mu$ l (to get this small volume an additional minute of centrifugation may be required).

## NOTES

**Estienne Swart** 30 Jan 2017

Recipe for complete sea water (1 L):

36 g Reef Crystals

1 ml Walne's solution

1 ml of 10  $\mu$ g/ml  $\text{FeSO}_4$

0.2 ml of 2 mg/ml thiamine (light sensitive; store at 4°C)

Add distilled water to 1 L

### Step 3.

Cells were transferred into two wells within a 96-well plates for transfection.

## ANNOTATIONS

**Rachele Cesaroni** 30 Jan 2017

One of the two wells was used for the negative control.

### Step 4.

FuGene HD Transfection Reagent was allowed to reach room temperature and mixed by inverting or vortexing briefly (if a precipitate is visible, briefly warm at 37°C and then let it reach room temperature).

#### **Step 5.**

90-98 µl of medium at room temperature (500 mM sorbitol, 0.5 mM Tris-HCl pH 7.0) and 4 µg of DNA in MilliQ H<sub>2</sub>O (0.2-1 µg/µL) were added to an Eppendorf tube and vortexed (after adding the DNA the final volume must be 100 µl).

#### **Step 6.**

6 µl of FuGene HD Transfection Reagent were added directly to the medium and mixed immediately.

#### **Step 7.**

FuGene HD Transfection Reagent and DNA mixture were incubated for 15 minutes at room temperature.

#### **Step 8.**

10 µl of the FuGene HD Transfection reagent and DNA mixture were added to the wells containing 50 µl of cells in medium (everything was mixed by pipetting).

#### **Step 9.**

One hour after addition of FuGene/DNA complexes, 50 µl of artificial sea water were added to the cells.

#### **Step 10.**

After another hour other 50 µl of artificial sea water were added to the cells.

#### **Step 11.**

An hour later cells were harvested (400 rcf for 3 minutes) and washed twice with artificial sea water (400 rcf for 3 min each time).

### **■ ANNOTATIONS**

**Rachele Cesaroni** 30 Jan 2017

Cells were resuspended in 400 µl of the artificial sea water.

#### **Step 12.**

Cells were transferred into glass depression wells for subsequent monitoring.

#### **Step 13.**

Cells were incubated at least for three days at 24°C, then examined by fluorescence microscopy to determine expression of the transformed construct.