



Oct 16, 2018

Working

Manual dissection of the *Schistosoma mansoni* head and back end for transcriptomic analysis [↗](#)

PLOS One

Alexandra Ley¹¹University Halle-Wittenberg

dx.doi.org/10.17504/protocols.io.taneide

Alexandra Ley

ABSTRACT

Total genomic DNA was extracted from leaf tissue to amplify and sequence individually five specific target loci from nuclear and chloroplast genome.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0205354>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Ley AC, Nissen J, Wölk A, Röser M (2018) Glacial refugia and speciation in a group of wind-pollinated and -dispersed, endemic Alpine species of *Helictotrichon* (Poaceae). PLoS ONE 13(10): e0205354. doi: [10.1371/journal.pone.0205354](https://doi.org/10.1371/journal.pone.0205354)

PROTOCOL STATUS

Working

1 Total genomic DNA was extracted from leaf tissue using the NucleoSpin®Plant-Kit (Macherey-Nagel, Düren, Germany) following the manufacturer's instructions.

2 Three chloroplast and one nuclear marker were amplified.

DNA-Region	Primer name	Sequence (5' → 3')	Reference
Chloroplast:			
<i>rps16</i>	<i>rpsF</i>	GTG GTA GAA AGC AAC GTG CGA CTT	[45]
	<i>rpsR2</i>	T CG GGA TCG AAC ATC AAT TGC AAT	[45]
<i>rpL32-tmL</i> (UAG)	<i>rpL32-F</i>	CAG TTC CAA AAA AAC GTA CTT C	[46]
	<i>tmL</i> (UAG)	CTG CTT CCT AAG AGC AGC GT	[46]
<i>ycf3ln1</i>	<i>ycf3ln1-F</i>	TGA CAG ATC ACG GCC ATA TT	[47]
	<i>ycf3ln1-R</i>	TTA YAG AGA TGG TGC GAT TT	[47]
<i>ycf3ln2</i>	<i>ycf3ln2-F</i>	GCT YGT TTC CAA TAC TCA GCA	[47]
	<i>ycf3ln2-R</i>	ATG GCC GTG ATC TGT CAT TA	[47]
Nuclear:			

<i>At103</i>	<i>At103-F</i>	CTT CAA GCC MAA GTT CAT CTT CTA	[48]
	<i>At103-R</i>	TTG GCA ATC ATT GAG GTA CAT NGT MAC ATA	[48]

- 3 Amplification of the target loci were was conducted in a Mastercycler (Eppendorf, Hamburg, Germany). Each 20 µl volume contains 2.00 µl 10× PCR buffer (without MgCl₂), 1.00 µl DMSO, 0.20 µl 100 mM dNTPs, 1.4 µl 50 mM MgCl₂, 0.20 µl 0.05 mM each forward and reverse primers, 0.20 µl Taq DNA polymerase (BioTherm DNA polymerase 5 u/µl from GeneCraft, Lüdinghausen, Germany), 1 to 2 µl genomic DNA extract and filled up with H₂O. Amplification cycles were as follows: one cycle of 3 min at 94°C, 35 cycles of 30 s at 94°C, 30 s at 50°C, 1.5 min at 72°C with a final extension period of 10 min at 72°C.
- 4 PCR products were sequenced by LGC Genomics (Berlin, Germany).
- 5 Sequences were edited and then aligned manually (after preliminary automatic alignment in MUSCLE [23]) with the program Geneious 4.8.3 (<http://www.geneious.com> [24]).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited