



Feb 26, 2019

Working

## Thawing and Seeding Frozen Cells

Kenneth Schackart<sup>1</sup><sup>1</sup>University of Arizona[dx.doi.org/10.17504/protocols.io.xkmfku6](https://doi.org/10.17504/protocols.io.xkmfku6) Kenneth Schackart 

### ABSTRACT

How to thaw cells from the liquid nitrogen storage and seed into a tissue culture flask

### PROTOCOL STATUS

#### Working

We use this protocol in our group and it is working

### GUIDELINES

Wear gloves. Spray anything entering or exiting the biosafety cabinet with 70% Ethanol.

### MATERIALS TEXT

- (1) T-75 or T-25 flask per frozen cell vial (or more if plating at a lower density)
- (1) 15 mL centrifuge tube per frozen cell vial
- (1) 10 mL serological pipette tip per T-75 flask or (1) 5 mL serological pipette tip per T-25 flask
- Warmed cell culture media
- 1000 µL filter pipette tips

### BEFORE STARTING

UV appropriate number of flasks, 15 mL centrifuge tubes, waste beaker, and serological pipette tips.  
Warm appropriate cell culture media.

#### Thaw Cells

- 1 Thaw cells by suspending in  37 °C water bath until completely thawed, but no longer than necessary

#### Transfer cell suspension

- 2 Within biosafety cabinet, transfer cell suspension to 15 mL centrifuge tube using 1000 µL pipette.

#### Dilute freezing medium

- 3 Add  1 ml warmed cell culture medium to cell suspension *dropwise*.



Adding the initial cell culture medium slowly helps prevent cell death caused by a rapid change in osmotic pressure.

- 4 Add an additional  3 ml warmed cell culture medium to cell suspension slowly.

#### Centrifuge cell suspension

#### Centrifuge Cell Suspension

- 5 Centrifuge the cell suspension at 1.5 kRPM for  00:03:00 .

#### Resuspend Cells

- 6 Remove bulk of supernatant with serological pipette, then remove remainder with 1000 µL pipette.



For small cell pellets, you are better off leaving a small amount of media than disturbing the cell pellet.

- 7 Add 1 mL warmed cell culture media to cell pellet.




Allowing the cell pellet rest in media for about 2 minutes will help with resuspension.

- 8 Gently pipette mix the cell pellet into the solution.

- 9 Add an additional 7 mL warmed cell culture media if using a T-75 flask.  
Add an additional 3 mL warmed cell culture media if using a T-25 flask.

#### Seed Cells

- 10 If using a T-75 flask, first add  2 ml warmed cell culture media to the flask.  
Using a serological pipette, transfer the cell suspension to the tissue culture flask.

#### Label Flask

- 11 Label the flask with:
- Cell line
  - Passage number
  - Date
  - Your initials

#### Incubate

- 12 Transfer flask to CO<sub>2</sub> incubator.

#### Documentation

- 13 Don't forget to remove the vial you used from the frozen storage inventory.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited