



Apr 16, 2019

Working

ATP assay [↗](#)

PLOS One

Kiichi Hirota¹, Yoshiyuki Matsuo¹¹Kansai Medical University[dx.doi.org/10.17504/protocols.io.v7ke9kw](https://doi.org/10.17504/protocols.io.v7ke9kw)

Yoshiyuki Matsuo

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215072>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Sumi C, Matsuo Y, Kusunoki M, Shoji T, Uba T, Iwai T, Bono H, Hirota K (2019) Cancerous phenotypes associated with hypoxia-inducible factors are not influenced by the volatile anesthetic isoflurane in renal cell carcinoma. PLoS ONE 14(4): e0215072. doi: [10.1371/journal.pone.0215072](https://doi.org/10.1371/journal.pone.0215072)

MATERIALS

NAME

CATALOG

VENDOR

CellTiter-Glo(R) Luminescent Cell Viability, 10ml

G7570

[Promega](#)

BEFORE STARTING

The CellTiter-Glo® luminescent cell viability assay kit (Promega) was used to evaluate the intracellular ATP content.

- 1 Cells were seeded in 96-well plates and allowed to grow for indicated time periods.
- 2 CellTiter-Glo reagent (50 µl) was then added directly into each well and incubated for 10 min prior to reading the plate using an EnSpire® Multimode Plate Reader (PerkinElmer, Waltham, MA, USA).

This detected the luminescence generated by the luciferase-catalyzed reaction between luciferin and ATP.
- 3 The ATP content was then calculated by comparing the luminescence levels of cells with those of control samples, with the latter defined as 100%.

Assays were performed in triplicate at least twice.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited