OPEN ACCESS



DNA extraction for vertebrate tissues using Phenol:Chloroform:Isoamylol

Chang Li

Abstract

This protocol is used to clarity the process of total DNA extration for L. maculatus genome.

Citation: Chang Li DNA extraction for vertebrate tissues using Phenol:Chloroform:Isoamylol. protocols.io

dx.doi.org/10.17504/protocols.io.ssyeefw

Published: 23 Aug 2018

Protocol

Sample preparation

Step 1.

- 1) Collect and cut muscle tissue of spotted sea bass.
- 2) The tissue was defrosted and washed with normal saline. Cut about 0.5g tissue.



REAGENTS

normal saline by <u>Sinopharm Chemical Reagent Co.</u>

Tissue lysis

Step 2.

Add STE, SDS (10%), protease K (20mg/ml) to mix. Heat 4-6h at 56° C, and shake in centrifuge at 500 rmp.

▼ TEMPERATURE

56 °C:



✓ STE by Contributed by users

SDS by Ambion

protease K by New England Biolabs

O DURATION

06:00:00 : heat

Phase separation

Step 3.

1) Add the same volume of Phenol : Chloroform : Isoamylol (25: 24: 1) and mix by gentle inversion for 3-5 min at room temperature.

- 2) Centrifuge at 13600rpm for 10min at 4° C and then transfer the supermatant aqueous phase in a new tube.
- 3) Add the same volume of Chloroform : Isoamylol (24 : 1) and mix by gently inverting for 3-5min at room temperature.
- 4) Centrifuge at 13600rpm for 10min at 4°C and then transfer the supermatant aqueous phase in a new tube.

```
    TEMPERATURE
4 °C: 2)
    TEMPERATURE
4 °C: 4)

REAGENTS
Phenol: Chloroform: Isoamylol (25: 24: 1) by Sinopharm Chemical Reagent Co.
Chloroform: Isoamylol (24: 1) by Sinopharm Chemical Reagent Co.
DURATION
00:05:00: 1)
DURATION
00:10:00: 2)
DURATION
00:05:00: 3)
DURATION
00:05:00: 3)
DURATION
```

00:10:00 : 4) DNA precipitation

Step 4.

- 1) Add equal volume of Isopropanol, and then mix by slowy inverting.
- 2) Incubate the tube at -20°C for overnight.
- 3) Centrifuge at 13600rpm for 10min at 4°C and discard the supermatant.

```
I TEMPERATURE
-20 °C: 2)
I TEMPERATURE
4 °C: 3)
I REAGENTS
I Isopropanol by Sinopharm Chemical Reagent Co.
UDURATION
00:10:00: 3)
```

DNA washing

Step 5.

- 1) Add 1ml of 70% of ethanol, flick the centriguge tube to resuspend the pellet, wash the DNA pellet for 3-5min.
- 2) Centrifuge at 13000rpm for 5min at 4°C, discard the ethanol.

- 3) Add 1ml of 70% ethanol and mix gently. Centrifuge at 13000rpm for 5min at 4°C, discard the ethanol.
- 4) Centrifuge at 13000rpm for 30s, discard the liquid. Air dry the pellet for 5-10min at room temperature.

```
■ TEMPERATURE
4 °C:2)
▮ TEMPERATURE
4 °C:3)
REAGENTS
5 70% of ethanol by Sinopharm Chemical Reagent Co.
O DURATION
00:05:00:1)
O DURATION
00:05:00:2)
O DURATION
00:05:00:3)
O DURATION
00:00:30:4)
O DURATION
00:10:00:4)
```

Step 6.

Dissolve DNA

Add RNase and 100ul of TE buffer (TE : RNase = 100:1), and then dissolve the pellet at 37° C for 30min in the Thermomixer.

