

# Full Contact Microbiology (a.k.a Diatom Transformation via Bacterial Conjugation) Version 2

Jeric Harper

## Abstract

**Citation:** Jeric Harper Full Contact Microbiology (a.k.a Diatom Transformation via Bacterial Conjugation). **protocols.io**  
<https://www.protocols.io/view/full-contact-microbiology-a-k-a-diatom-transformat-fzrbp56>

**Published:** 18 Oct 2016

## Before start

- Make sure you are using fresh *E. coli* cells streaked for isolation on LB + antibiotics no more than 1 week from -80°C cryostock.

## Protocol

### Diatom cultivation (liquid)

#### Step 1.

**Grow** the diatom culture to mid-log phase ( $\approx 8.0E^6$  cells/ml for *Phaeodactylum tricornutum* grown on F/2 media or  $5.0 E^7$  cell/ml when grown on BG-11). We have investigated transformation efficiency throughout the growth curve and found this to be the sweetspot.

#### 📌 NOTES

**Jeric Harper** 04 Oct 2016

#### Cultivation Conditions

- Cultivation in 0.2  $\mu$ m-filtered L1 or BG-11 media prepared using 32ppt seawater collected from the Gulf of Mexico.
- 80 - 100  $\mu$ mol photons  $m^{-2} s^{-1}$  provided by Philips Daylightt Delux 40 watt T12 flourscent tubes.
- Semi-continuous operation in 1L bubble column bioreactors (500 - 800 ml working volume).
- 0.2  $\mu$ m filtered air containing 1%  $CO_2$  provides aeration/agitation at 0.2 vvm.

#### ■ ANNOTATIONS

**Jeric Harper** 04 Oct 2016

- [The original protocol](#) plated the culture (250uL of 1.0E8 cells/ml) on 1/2 strength L1, 1% agar plates for 4 days prior to transformation.

- We have been investigating liquid cultivation because many diatom species do not survive on agar.

## Prepare *E. coli* starter culture

### Step 2.

**Grow** 1 mL of *E. coli* culture containing both the mobility plasmid (Pta-MOB) and carrier plasmid, overnight (16-20 hrs) in LB+antibiotics, for each planned transformation. (We grow them at 37°C at 270 rpm in a shaking incubator.)

#### DURATION

16:00:00

#### ANNOTATIONS

**Jernej Turnsek** 19 Oct 2016

In my experience it has proven beneficial to perform conjugations with multiple donors carrying THE SAME construct. 2/3 worked beautifully, one failed completely.

## Outgrow *E. coli*

### Step 3.

On the day of transformation, use the overnight culture to inoculate 50 mL of fresh LB+antibiotic, 1:50 dilution, for each planned transformation.

- Grow to an OD<sub>600</sub> of 0.8 - 1.0 (37°C with 270 rpm shaking).
- This takes about 3-4 hours.

#### DURATION

03:00:00

#### ANNOTATIONS

**Jeric Harper** 06 Oct 2016

The authors of [The original protocol](#) reported that the OD<sub>600</sub> range is flexible. Transformation success has been seen within OD<sub>600</sub> range of 0.4 to 1.2.

## *P. tricornutum* cell concentration

### Step 4.

During the 3-4 hours the *E. coli* culture is growing, measure the *Phaeodactylum tricornutum* cell concentration with a FlowCam or haemocytometer to calculate the required volume needed to collect 2.5E<sup>8</sup> cells for each transformation.

#### DURATION

00:05:00

## Concentrate the diatom and *E. coli* cultures

## Step 5.

For each transformation, centrifuge 50 mL of *E. coli* culture and the required *Phaeodactylum tricornutum* volume at 4000 x g for 10 minutes at 4°C.

- Resuspend *E. coli* pellet in 500 µL of SOC medium.
- Resuspend *P. tricornutum* pellet in 500 µL of L1 medium.

**Note:** The diatom and *E. coli* cultures should be **centrifuged at around the same time** to minimize the amount of time they spend concentrated.

### DURATION

00:10:00

### NOTES

**Jeric Harper** 04 Oct 2016

[The original protocol](#) scraped the agar plates the culture was initiated on using 500uL F/2, then adjusted the volume to attain  $5.0 \times 10^8$  cells/ml. We have found no difference so far in transformation efficiency between liquid- and plate-initiated cultures.

## Conjugation

### Step 6.

In a 1.5 mL tube mix 200 µl of *E. coli* cells with 200 µl of *Phaeodactylum tricornutum* cells.

Negative control: In a 1.5 mL tube mix 200 µl of SOC medium with 200 µl of *Phaeodactylum tricornutum* cells.

**Note:** Incubate and treat the negative control plates identically to conjugation plates.

### ANNOTATIONS

**Jeric Harper** 04 Oct 2016

- The authors of [The original protocol](#) suggested spreading the mixture near, but not touching, the edges of the plate, as this mixture will be scraped and replated in 2 days.
- N concentration seems to be important. We have found greater transformation success in using 0.5x BG-11 medium made with seawater instead of 0.5x L1 medium.

- 5% LB is a sufficient concentration to support the *E. coli* culture without encouraging growth. Increasing the LB concentration results in increased algal death. We have also buffered the media but that did not increase *P. tricornutum* survivorship.

## Conjugation

### Step 7.

Spread the mixture (400 µL) on [Conjugation Plates](#). (0.5x BG-11 with 5% LB and 1% agar).

#### 📌 NOTES

**Jeric Harper** 06 Oct 2016

- The authors of [The original protocol](#) suggested spreading the mixture near, but not touching, the edges of the plate, as this mixture will be scraped and replated in 2 days.
- N concentration seems to be important. We have found greater transformation success in using 0.5x BG-11 medium made with seawater instead of 0.5x L1 medium.
- 5% LB is a sufficient concentration to support the *E. coli* culture without encouraging growth. Increasing the LB concentration results in increased algal death. We have also buffered the media but that did not increase *P. tricornutum* survivorship.

## Conjugation

### Step 8.

Incubate plates for 90 minutes at 30°C in the dark.

#### 🕒 DURATION

01:30:00

#### 📌 NOTES

**Jeric Harper** 06 Oct 2016

*P. tricornutum* can survive transformation temperatures up to 32°C. At 34°C and above survivorship is severely hampered.

## Conjugation

### Step 9.

Move plates to light incubator (18°C and 100 µmol photons m<sup>-2</sup> s<sup>-1</sup>) for 2 days.

## ⌚ DURATION

48:00:00

## 📌 NOTES

**Jeric Harper** 06 Oct 2016

- According to [Diner et. al](#), the conjugation occurs during the 2-day incubation. Transformation efficiency increased as incubation time increased

## Selection

### Step 10.

Collect cells by adding 1 mL of L1 medium. Use a cell scraper to concentrate cells and medium to one side of the plate. Transfer resuspended cells to a 1.5 ml microcentrifuge tube with a P1000 pipette and filter tips.



## ■ ANNOTATIONS

**Jeric Harper** 04 Oct 2016

It is important to use filter tips as small pieces of agar that are accidentally scraped up can cause the pipette to cavitate resulting in contamination.

We will often first add 500uL to collect the majority of the cells from the plate into the 1.5mL tube, then repeat with an additional 500ul to collect the remainder.

## Selection

### Step 11.

Spread 200 µl of the cell suspension on a [Selection Plate](#).

## 📌 NOTES

**Jeric Harper** 06 Oct 2016

- Dilution of the selection volume (using L1 as makeup volume) can help facilitate colony enumeration if 200µl results in too many colonies.
- To maximize the number of colonies the entire volume of resuspended cells can be plated to multiple plates.
- Collect 1 µL of the cell suspension and dilute 1:1000 for FlowCam cell count.

## Selection


### Step 12.

Incubate at 18°C and 100 µmol photons m<sup>-2</sup> s<sup>-1</sup> until colonies appear.

## Colony identification

### Step 13.

After a minimum of 8-12 days, untransformed *Phaeodactylum tricornutum* cells die off, and colonies of transformed cells begin to appear – in some cases, this can take 3-4 weeks.




Alternatively selection can be done in liquid BG-11 [Selection media](#) using eGFP as a reporter and sorted using FACS.

un-transformed Pt



eGFP expression



#### ■ ANNOTATIONS

**Jeric Harper** 06 Oct 2016

- For plate selection, use [ImageJ protocol](#) for colony enumeration.

- Calculate transformation efficiency using:

Efficiency = (number of colonies on plate) ÷ (selection volume cell density (cell/mL) x volume (mL) put on selection plate)

**Note:** Be sure to include a dilution factor in the calculation if it was used in plating or counting!

#### Step 14.

This protocol was modified from [the original procedure](#) and correspondence with the authors.