

Propidium Iodide (PI) Staining Method

Dr. Steven Wilhelm

Abstract

Please contact Dr. Steven Wilhelm (wilhelm@utk.edu) for additional information regarding this protocol.

Modified from J Mol Biol 13,269 (1965)

Citation: Dr. Steven Wilhelm Propidium Iodide (PI) Staining Method. **protocols.io**

dx.doi.org/10.17504/protocols.io.ibxcapn

Published: 13 Jun 2017

Protocol

Sample Preparation

Step 1.

Make a sample by adding 25% glycerol to cells



REAGENTS

Glycerol MRGE-4002 by growcells.com



NOTES

Alyssa Alsante 06 Jun 2017

Samples can be stored at -20°C for 30 days.

Staining Solutions Preparation

Step 2.

Add DAPI at 50 µg/mL in dH₂O

Staining Solutions Preparation

Step 3.

Add PI at 500 µg/mL in dH₂O

Staining Solutions Preparation

Step 4.

Filter sterilize with 0.2 µm Millipore filter and store in the dark at 4°C

Slide Preparation

Step 5.

Dilute cells appropriately according to your sample

Step 6.

Add DAPI at a final concentration of 5 µg/mL

Step 7.

Add PI at a final concentration of 3 µg/mL

Step 8.

Incubate in the dark for 30 min

 DURATION

00:30:00

Step 9.

Filter sample onto a 0.2 µm black polycarbonate filter under low vacuum (<15 mmHg)

Step 10.

Mount the filter onto a glass slide with low fluorescence immersion oil

Step 11.

View the sample under the Leica DM4000-6000 epifluorescent microscope under the Texas Red N3 filter at green excitation = 546 nm and emission = 600 nm