

# Molluscum contagiosum real-time PCR

Judy Northill, Russell Simmons

## Abstract

A real-time PCR method for the detection of Molluscum contagiosum virus from human samples. This assay targets the MC021L gene and detects both subtype 1 and 2 of the virus.

**Citation:** Judy Northill, Russell Simmons Molluscum contagiosum real-time PCR. **protocols.io**

dx.doi.org/10.17504/protocols.io.kksuwe

**Published:** 06 Nov 2017

## Before start

- If using a different brand or model of real-time thermocycler, check the concentration of ROX is adequate.
- Method assumes the user is familiar with the thermocycler and software used to run the protocol.

## Protocol

### Oligonucleotide sequences

#### Step 1.

Name	5'-3'
MCVp43kF (forward primer)	GCTCACGTACGACTGCTTYGAC
MCVp43kR (reverse primer)	CGTGGAGCGCAGATTGC
MCVp43kP (probe)	6FAM-CGCTCATCTCGCAGAC-MGB

### Reaction set-up

#### Step 2.

Assay has been used on both a Rotor-Gene 6000 / Rotor-Gene Q 5-plex using 100-place rotor discs and an ABI 7500 using 96-well plates.

Total reaction volume is 20µL and is suitable for both formats.

Prepare sufficient for number of reaction plus a 'dead volume' usually 2 extra. Adjust as necessary if using a robotic dispenser.

Reagent	Vol ( $\mu$ L) x1	Final reaction concentration
Nuclease-free water	4.91	
MCVp43kF 200pmol/ $\mu$ L	0.03	300nM
MCVp43kR 200pmol/ $\mu$ L	0.03	300nM
MCVp43kP 100pmol/ $\mu$ L	0.03	150nM
<sup>1</sup> TaqMan™ Fast Universal PCR Master Mix (2X)	10	1X
<b>Template</b>	<b>5</b>	

<sup>1</sup>Thermofisher product [4352042](#)

Dispense 15 $\mu$ L to each reaction well.

Add 5 $\mu$ L of template, extracted DNA, controls or NTC (nuclease-free water).

Total reaction volume is 20 $\mu$ L



## REAGENTS

TaqMan™ Fast Universal PCR Master Mix (2X) [4352042](#) by [Applied Biosystems](#)

## Amplification

### Step 3.

The assay has been optimised and validated for the ABI 7500. It also is used on the Rotor-Gene 6000 and Rotor-Gene Q thermocyclers.

## PCR

50°C 5min

95°C 2min

95°C 3s 40X

60°C 30s\*

\*Florescence acquisition step

## Result analysis

### Step 4.

The threshold should be placed in the exponential range above any background noise within the assay.

A positive result is one where the  $C_T$  is  $<40$  and produces a sigmoidal curve.

NTC should not produce a curve and should be greater than  $40C_T$ .