

# Iron Chloride Precipitation of Viruses from Seawater

Seth John, Bonnie Poulos, Christine Schirmer

## Abstract

This protocol describes a technique to recover viruses from natural waters using iron- based flocculation and large-pore-size filtration, followed by resuspension of virus- containing precipitates in a pH 6.5 buffer. This Fe-based virus flocculation, filtration and resuspension method (FFR) is efficient (> 90% recovery), reliable, inexpensive and adaptable to many aspects of marine viral ecology and genomics research. Recovered viruses are amenable to gene sequencing, and a variable proportion of phages, depending upon the phage, retain their infectivity when recovered if using oxalic acid in the resuspension buffer. Particles lose infectivity if resuspended with ascorbic acid in the buffer.

**Citation:** Seth John, Bonnie Poulos, Christine Schirmer Iron Chloride Precipitation of Viruses from Seawater. **protocols.io**

dx.doi.org/10.17504/protocols.io.c2wyfd

**Published:** 20 Jul 2015

## Guidelines

**PCR inhibition:** Note, that due to the likely issue with the EDTA in the resuspension buffer interacting with the MgCl in a PCR reaction, direct PCR is ineffective. To circumvent this, extract DNA using Wizard columns from Promega (see [Tucson Marine Phage Lab's DNA Extraction of Viruses Using Wizard Columns protocol](#)).

### Materials:

- Seawater
- Tris-base (FW=121.14)
- Na<sub>2</sub>-EDTA dihydrate (FW=372.24)
- MgCl<sub>2</sub> hexahydrate (FW=203.3)
- Ascorbic acid (FW=176.12) and/or Oxalic acid (FW=126.07)
- NaOH, 5N
- HCl
- FeCl<sub>3</sub>·6H<sub>2</sub>O (FW=270.3)
- MilliQ water
- pH paper
- Whatman glass fiber filter GF/A (cat. 1820-150; 1.6µm retention; 150mm diameter)
- Millipore Express PLUS filter (cat. GPWP14250; 142mm, 0.22µm)
- Pall Supor-800 filter (cat. 60114; 142mm, 0.8µm)
- GE Polycarbonate Membrane filter (cat. K10CP14220; 142mm, 1.0µm), Available through Midland Scientific Inc. (cat. MAINE 1216611)
- pH meter
- 142mm filtering towers (Need 2 for efficiency)
- Peristaltic pump with a pressure gauge

## Notes:

**Alternative Pre-filtering Method:** An alternative method for preparing seawater is to pre-filter the seawater through a Whatman GF/D glass fiber filter (cat. 1823-150; 2.7µm retention; 142mm diameter) to remove large particles. The filtrate is filtered again using a Millipore Steripak GP10 or GP20 (cat. SPGPM10RJ or SPGPM20RJ; 0.22µm). The GP 10 works for <10L volumes and GP20 works for <20L volumes. Larger capsule filters will easily filter 100-200 liters: Pall Acropak 200, 500, 1000 and 1500. The cellular fraction is retained on the filter while the virus fraction is in the filtrate.

Filtration should be done at a maximum pressure of 15 psi. A peristaltic pump applies the pressure and a stainless steel Millipore filter rig (PN YY3014236) holds filters. For less expense, use an acrylic or polycarbonate 142 mm filter holder (Geotech Environmental Equipment, acrylic or PC) and pressurize carboy with an air pump. Most lab vacuum pumps have an air-pressure outlet on the other side that can easily be adjusted to provide 15 psi pressure. Nalgene sells caps for their 20L carboys with fittings in the top that can be used for pressurization.

**FeCl<sub>3</sub> Formulation:** The solution of ferric chloride hexahydrate is calculated based on the amount of iron, not on the amount of the salt. The stock solution has 10g Fe per liter. For precipitation, the final optimal concentration of Fe in seawater (SW) is 1mg Fe per liter SW which is equal to 2.9mg FeCl<sub>3</sub> per liter SW or to 4.83mg FeCl<sub>3</sub>·6H<sub>2</sub>O per liter SW.

**Equipment Care:** Carboys and collection bottles should be washed with 1N HCl and rinsed with MilliQ water before use. Checking tubing before use and replace if worn. Always rinse stainless steel filter rigs with MilliQ water after use and air dry to prevent corrosion.

**Original Author(s):** Seth John, Bonnie Poulos

**Revision Author(s):** Bonnie Poulos, Christine Schirmer

**Citation:** John, S.G., Mendez, C.B., Deng, L., Poulos, B., Kauffman, A.K.M., Kern, S., Brum, J., Polz, M.F., Boyle, E.A., & Sullivan, M.B. (2011). A simple and efficient method for concentration of ocean viruses by chemical flocculation. *Environ Microbiol Rep.* 3(2), 195-202.  
doi:10.1111/j.1758-2229.2010.00208.x.

## Materials

- ✓ Millipore Express PLUS filter 0.22 µm, 142 mm [GPWP14250](#) by Contributed by users
- ✓ Pall Supor-800 filter [60114](#) by Contributed by users
- ✓ GE Polycarbonate Membrane filter [MAINE 1216611](#) by Contributed by users

## Protocol

### Virus Precipitation

#### Step 1.

Assemble two 142mm filtering towers and attach the filter apparatus to a peristaltic pump with a pressure gauge (maximum pressure = 15 psi)

#### NOTES

**Lauren Chittick** 29 Jun 2015

See Guidelines note for an alternative pre-filtration method.

#### ■ ANNOTATIONS

**Lauren Chittick** 17 Feb 2016

For smaller volume testing, we used 1.0-1.5L filtrate per 47mm tower.

### Virus Precipitation

#### Step 2.

Pre-filter 20L of seawater using a 150mm Whatman GF/A filter.

#### 🔗 NOTES

**Lauren Chittick** 30 Jun 2015

Wear gloves and use forceps for handling all filters.

### Virus Precipitation

#### Step 3.

Filter again to remove bacteria using a 0.22µm, 142mm Millipore Express Plus filter.

#### 🧴 REAGENTS

✓ Millipore Express PLUS filter 0.22 µm, 142 mm [GPWP14250](#) by Contributed by users

#### 🔗 NOTES

**Lauren Chittick** 23 Jun 2015

See note in guidelines for alternative pre-filtration method.

**Lauren Chittick** 30 Jun 2015

Wear gloves and use forceps for handling all filters.

### Virus Precipitation

#### Step 4.

Treat the virus fraction (0.22µm filtrate) with FeCl<sub>3</sub> to precipitate the viruses by adding 1ml of 10g/L Fe Stock Solution to each 20L of filtrate. Shake vigorously for 1 minute.

#### 📄 AMOUNT

1 ml Additional info:

#### 🕒 DURATION

00:01:00

#### 📋 PROTOCOL

#### . [10g/L Fe Stock Solution](#)

CONTACT: [Bonnie Poulos](#)

#### 🔗 NOTES

**Lauren Chittick** 13 Jul 2015

Solution is acidic and should be handled with care.

Discard solution if a cloudy precipitate forms. Do not dilute solution as iron hydroxide precipitate will form quickly.

For use on board ship, it is best to pre-weigh iron chloride into 50ml centrifuge tubes and add water as needed during the cruise: pre-weigh 0.966g FeCl<sub>3</sub>·6H<sub>2</sub>O into several 50ml tubes. Add 20ml MQ-H<sub>2</sub>O when ready to precipitate viruses in seawater.

#### Step 4.1.

Measure out 4.83g  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ .

 [AMOUNT](#)

5 g Additional info:

 [REAGENTS](#)

✓ Iron(III) chloride hexahydrate [236489](#) by Contributed by users

#### Step 4.2.

Dissolve  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  in 100ml MilliQ water.

 [AMOUNT](#)

100 ml Additional info:

#### Step 4.3.

Store at room temperature or 4°C.

### Virus Precipitation

#### Step 5.

Add an additional 1ml of 10g/L Fe stock solution to each 20L of filtrate (for a total of 2ml Fe Stock Solution per 20L filtrate). Shake vigorously for 1 minute.

 [AMOUNT](#)

1 ml Additional info:

 [DURATION](#)

00:01:00

 [PROTOCOL](#)

#### . [10g/L Fe Stock Solution](#)

CONTACT: [Bonnie Poulos](#)

#### [NOTES](#)

**Lauren Chittick** 30 Jun 2015

For use on board ship, it is best to pre-weigh iron chloride into 50ml centrifuge tubes and add water as needed during the cruise: pre-weigh 0.966g  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  into several 50ml tubes. Add 20ml MQ-H<sub>2</sub>O when ready to precipitate viruses in seawater.

**Lauren Chittick** 30 Jun 2015

Discard solution if a cloudy precipitate forms. Do not dilute solution as iron hydroxide precipitate will form quickly.

**Lauren Chittick** 30 Jun 2015

Solution is acidic and should be handled with care.

#### Step 5.1.

Measure out 4.83g  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ .

 [AMOUNT](#)

5 g Additional info:

 [REAGENTS](#)

✓ Iron(III) chloride hexahydrate [236489](#) by Contributed by users

#### Step 5.2.

Dissolve  $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$  in 100ml MilliQ water.

 [AMOUNT](#)

100 ml Additional info:

#### Step 5.3.

Store at room temperature or 4°C.

## Virus Precipitation

### Step 6.

Let the  $\text{FeCl}_3$  treated filtrate sit for 1 hour at room temperature.

#### DURATION

01:00:00

## Virus Precipitation

### Step 7.

Filter the  $\text{FeCl}_3$  treated filtrate using a 1.0 $\mu\text{m}$ , 142mm, polycarbonate (PC) membrane filter on top of a 0.8 $\mu\text{m}$ , 142mm, Supor support filter and attaching the filter apparatus to a peristaltic pump with a pressure gauge (maximum pressure = 15 psi).

#### REAGENTS

✓ GE Polycarbonate Membrane filter [MAINE 1216611](#) by Contributed by users

#### NOTES

**Lauren Chittick** 30 Jun 2015

The  $\text{FeCl}_3$  precipitate will be captured on the PC filter; the Supor is just for support in the 142mm stainless steel filtration apparatus. Although the 20L can be collected on a single 142mm PC membrane, it is faster to change the PC membrane one or two times during the filtration; normally 3 PC membranes and 1 Supor membrane are used per 20L seawater.

## Virus Precipitation

### Step 8.

Place all of the polycarbonate filters from the 20L into a 50ml centrifuge tube being careful not to scrape off any of the  $\text{FeCl}_3$  on the edge of the tube (having precipitate facing out aids in dissolving the precipitate).

## Virus Precipitation

### Step 9.

Discard the Supor support filter.

#### NOTES

**Lauren Chittick** 30 Jun 2015

The filters can be stored in the 50ml tubes at 4°C until ready to resuspend the precipitated viruses. Be sure the caps are on securely so that the filters do not dry out.

## Resuspension

### Step 10.

Prepare fresh 0.1M EDTA-0.2M  $\text{MgCl}_2$ -0.2M Ascorbate Buffer, pH6.5.

#### PROTOCOL

. [0.1M EDTA-0.2M  \$\text{MgCl}\_2\$ -0.2M Ascorbate Buffer](#)

CONTACT: [Bonnie Poulos](#)

#### ANNOTATIONS

**Lauren Chittick** 30 Jun 2015

Fresh buffer with no precipitate and a pH of 6.0 are critical!

**Lauren Chittick** 20 Jun 2016

For an easier formulation, see the protocol for making 1M Citrate and Magnesium Resuspension

Buffer.

#### 1x Buffer

##### Step 10.1.

Dissolve 1.51g Tris-base in 80ml Milli Q water.

#### 1x Buffer

##### Step 10.2.

Dissolve 3.72g Na<sub>2</sub>-EDTA dihydrate into solution.

#### 🔗 NOTES

**Bonnie Poulos** 15 Jun 2015

pH will be ~10.0

#### 1x Buffer

##### Step 10.3.

Once EDTA is in solution, dissolve 4.07g MgCl<sub>2</sub>.

#### 🔗 NOTES

**Bonnie Poulos** 15 Jun 2015

pH will drop to ~8.0

#### 1x Buffer

##### Step 10.4.

Add 3ml of NaOH.

#### 🔗 NOTES

**Bonnie Poulos** 15 Jun 2015

This will drop the pH to ~4.5 and the solution will become cloudy which indicates that the EDTA is coming out of solution.

#### 1x Buffer

##### Step 10.5.

Dissolve the reductant (3.52g of ascorbic acid or 2.52g of oxalic acid).

#### 🔗 NOTES

**Bonnie Poulos** 15 Jun 2015

The pH will increase to ~8.3 and the solution will clear up.

#### 📌 ANNOTATIONS

**Uri Neri** 10 Apr 2018

Dear Bonnie,

We're trying to create the resuspension buffer for the VLPs-Iron precipitates resuspension step, with the reductant agent being oxalic acid (anhydrous, FW=90.04). net weight for the acid concentration was calculated to adjust for the anhydrosity and is ~1.808 [g]. Other than this minor change, I've made no modifications to the protocol, yet I am unable to create the buffer without the solution becoming 'murky' after the addition of 3ml NaOH(5N), and after some time without stirring, visible precipitates are formed (and the solution becomes clear).

We've followed the protocol to the letter, and continuously measured pH.

Your thoughts on the matter would be greatly appreciated.

With best regards,

Uri Neri

1x Buffer

#### Step 10.6.

Once the reductant is in solution, add the last 1ml of NaOH.

1x Buffer

#### Step 10.7.

Check the pH using pH paper (the buffer should be at pH 6.0 - 6.5)

#### 🔗 NOTES

**Bonnie Poulos** 23 Jun 2015

The solution may need some minor adjusting with NaOH or HCl to achieve a pH of 6.0.

**Bonnie Poulos** 23 Jun 2015

pH 6.0 is ideal for good recovery of viruses.

1x Buffer

#### Step 10.8.

Check the volume and add MilliQ water for a total volume of 100ml.

1x Buffer

#### Step 10.9.

Store the buffer in the dark (bottle wrapped in foil) and visually inspect prior to use. It should be clear without precipitates.

#### 🔗 NOTES

**Bonnie Poulos** 15 Jun 2015

At this point, 10-15ml of buffer can be sacrificed for a final pH check using a pH meter.

**Bonnie Poulos** 23 Jun 2015

The buffer will start to change color after about 24 hours. It is okay to use if slightly discolored, but do not use after about 36 hours (eventually the buffer will turn almost orange!).

### Resuspension

#### Step 11.

Add 1ml buffer for each 1mg Fe, which is 20ml of 1x buffer for the 20L seawater precipitate (10ml of a 2X concentrated buffer). Shake the tube vigorously.

### Resuspension

#### Step 12.

Place tubes on a rotator and rotate at 4°C overnight.

#### 🕒 DURATION

18:00:00

#### 🔗 NOTES

**Lauren Chittick** 30 Jun 2015

You will see the solution change colors a number of times – that is okay. However, if precipitates form, there is a possibility that viruses are trapped or adhered to the particles. Centrifuge gently to pellet the precipitate, and try using additional buffer to get the pellet back into solution. Although it may not go back into solution, the pellets will have been washed and may recover some viruses.

13. Leave stored at 4°C until ready to analyze or further purify the viral fraction.

### Resuspension

#### Step 13.

Leave stored at 4°C until ready to analyze or further purify the viral fraction.

#### Resuspension

##### Step 14.

To collect the resuspended precipitate, open the tube and collect liquid into a fresh tube.

#### Resuspension

##### Step 15.

Using a sterile pipet, applicator stick, or bleached and rinsed forceps pull an edge of the filter onto the tube edge so it is caught in the lid when secured.

#### Resuspension

##### Step 16.

Centrifuge at low speed (less than 1000 rpm) for 3-5 minutes to recover more buffer.



DURATION

00:05:00

#### Resuspension

##### Step 17.

If the filter has precipitate on it, add more buffer, incubate again for a couple hours, and collect liquid as stated above.

#### ■ ANNOTATIONS

**Lauren Chittick** 17 Feb 2016

Citrate buffers were a little more "stubborn" so I briefly vortexed and rinsed with the collected buffer (instead of adding new buffer).

## Warnings

Do NOT let the filters dry out before the precipitated viruses are resuspended! This will cause a loss in yield. - It is important that the filters remain moist until resuspension. The excess liquid on top of the filter can provide enough humidity if the tubes are well sealed with a tight-fitting cap and parafilm. One or two ml of sterile water or 0.2um filtered SW from the same source can be added if necessary to keep the filters moist. - Keep filters at 4°C until processed, and processing is best if done soon after collection (although we have successfully resuspended particles after several years if stored properly).