

Mar 28,  
2019

Working

Chromatographic profile of aqueous extract by HPLC [↗](#)

Version 3

PLOS One

Luis Wiliunfo Torres Tapia<sup>1</sup><sup>1</sup>Centro de Investigación Científica de Yucatán A. C.[dx.doi.org/10.17504/protocols.io.wy9ffz6](https://doi.org/10.17504/protocols.io.wy9ffz6) **Luis Wiliunfo Torres Tapia**  
Centro de Investigación Científica de Yucatán A. C. 

## ABSTRACT

This chromatographic detection method was developed in order to solve the most compounds from aqueous extracts leaves of dune and mangrove plant by HPLC.

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0213493>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Rodríguez-García CM, Ruiz-Ruiz JC, Peraza-Echeverría L, Peraza-Sánchez SR, Torres-Tapia LW, Pérez-Brito D, Tapia-Tussell R, Herrera-Chalé FG, Segura-Campos MR, Quijano-Ramayo A, Ramón-Sierra JM, Ortiz-Vázquez E (2019) Antioxidant, antihypertensive, anti-hyperglycemic, and antimicrobial activity of aqueous extracts from twelve native plants of the Yucatan coast. PLoS ONE 14(3): e0213493. doi: [10.1371/journal.pone.0213493](https://doi.org/10.1371/journal.pone.0213493)

## PROTOCOL STATUS

Working

## MATERIALS

NAME	CATALOG #	VENDOR
Trifluoroacetic acid for HPLC > 99.0%	302031-100ML	<a href="#">Sigma-aldrich</a>
Acetonitrile HPLC	9012-03	<a href="#">fisher</a>
Methanol HPLC	9093-03	<a href="#">Fisher Scientific</a>

## SAFETY WARNINGS

The trifluoroacetic acid is very high volatile, try to use a extraction cabinet to take it.

## BEFORE STARTING

Everything that goes into HPLC must be filtered first, through a 0.2 µm nylon filter and special glassware to remove particles that can get caught up on the column to avoid overpressure problems.

The potentiometer must be calibrated.

## Sample preparation

- 1 Each aqueou extract was added a 10% MeOH to improve solubility, and next it was filtered through 0.2 µm filter disc.

## Movil phase

- 2 Movil phase consist:
- Movil phase A contained ultrapure type 1 water ( Simplicity® Water Purification System, Millipore) adjusted to pH 2.5 with trifluoroacetic acid (TFA).
  - Movil phase B contained acetonitrile (ACN).
- Each movil phase was filtered through 0.2 µm filter and finally was degassed during 30 min. in ultrasonic equipment.

#### Equipment

- 3 - Profile was carried out on an Agilent Series 1290 Infinity HPLC System with vacuum degasser, quaternary pump, autosampler, thermostated column compartment and photodiode array detector (DAD).
- Data analysis was performed with Agilent HPLC EZChrom software.

#### Column

- 4 Separation of analytes was performed on an Grace Alltima C18 column (4.6 mm × 250 mm, 5 µm), which was maintained at 35 °C during the analysis.

#### Elution gradient

- 5 The gradient was programmed as follows: 0-3 min, 5% B; 3-43 min, 5-30% B; 43-73 min, 30-85 %B, 73-75 min 85-5 %B.

#### Elution flow

- 6 Set the flow rate to 1 mL/min.

#### DAD monitoring

- 7 Set the DAD monitoring to 254 nm and 350 nm.

#### Injection

- 8 Inject 20 µL of the sample.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited