



Transfection of constructs in diplomids to enhance the HR pathway, using RS-1 (3-(benzylamino) sulfonyl)-4-bromo-N-(4-bromophenyl) benzamide).

Version 3

Binnypreet Kaur¹, Drahomíra Faktorová², Julius Lukeš³

¹ [Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic], ² [Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic], ³ [Institute of Parasitology Czech Academy of Sciences Faculty of Sciences University of South Bohemia Czech Republic] [Canadian Institute for Advanced Research Toronto Canada]

dx.doi.org/10.17504/protocols.io.x45fqy6

Protist Research to Optimize Tools in Genetics (PROT-G)

Julius Lukes

 Binnypreet Kaur   

ABSTRACT

RS-1 (RAD51-stimulatory compound 1) is a stimulator of homologous recombination (HR) protein RAD51

PROTOCOL STATUS

In development

We are still developing and optimizing this protocol

- 1 Step 1: Count the cells and plan to 1 nucleofection with 5×10^7 cells for each construct.
- 2 Step 2: Harvest the cells by centrifugation at 1300xg for 5 min at room temperature in Swing Bucket Rotor.
- 3 Step 3: Resuspend the cell pellet in 100ul of AMAXA Human T- cell solution at 4C (from refrigerator combine 81.8ul of Human T-cell nucleofector solution + 18.2ul Supplement).
- 4 Step 4: Add 5-10ug of (PCR) (linearized DNA) along with 0.3 uM RS-1 enhancer into the cuvette (resuspend in 10ul of H2O).
- 5 Step 5: Put everything into the cuvette, close the cap and place in the electroporator, cuvette should only fit in one direction, but metal sides should face towards you.
- 6 Step 6: Press for the Program X-001 to electroporate and repeat the same as mentioned in transformation protocol.
- 7



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits

