

# Seydoux lab Cas9 preparation

Alexandre Paix, Andrew Folkmann, Dominique Rasoloson, Geraldine Seydoux

## Abstract

This protocol is for purification of Cas9::NLS<sub>SV40</sub>::His<sub>6</sub> from:

Paix A, Folkmann A, Rasoloson D, and Seydoux G (2015) [High Efficiency, Homology-Directed Genome Editing in Caenorhabditis elegans Using CRISPR/Cas9 Ribonucleoprotein Complexes](#). Genetics genetics.115.179382; Early online July 17, 2015. doi:10.1534/genetics.115.179382

Please see the [full manuscript](#) for detail.

**Citation:** Alexandre Paix, Andrew Folkmann, Dominique Rasoloson, Geraldine Seydoux Seydoux lab Cas9 preparation. **protocols.io**

dx.doi.org/10.17504/protocols.io.dn55g5

**Published:** 06 Oct 2015

## Guidelines

Protocol updates will be posted on the [Seydoux lab website](#).

Main protocol: [Direct delivery CRISPR-HDR editing protocol for C. elegans](#)

### Buffers:


Buffer A: 20mM Tris pH 8.0, 250 mM KCl, 20 mM imidazole, 10% glycerol, 1 mM TCEP

Buffer B: 20mM Tris pH 8.0, 800 mM KCl, 20 mM imidazole, 10% glycerol, 1mM TCEP

Buffer C: 20mM Hepes pH 8.0, 500 mM KCl, 250 mM imidazole, 10% glycerol

Buffer D: 20mM Hepes pH 8.0, 500 mM KCl, 20% glycerol

### Purified Cas9::NLS<sub>SV40</sub>::His<sub>6</sub> resolved by SDS-PAGE:



Recombinant Cas9::NLS<sub>SV40</sub>::His<sub>6</sub> was affinity purified using Ni<sup>2+</sup> agarose (lane 1). Pooled eluent was flowed over Q sepharose to remove contaminating DNA bound to Cas9 (lane 2). Samples were resolved by SDS-PAGE and visualized by coomassie staining.

**Cas9 activity assay:** We recommend testing your Cas9 preparation using [this method](#).

## Protocol

### Step 1.

Transform DE3 GOLD (Agilent, #230132) cells with nm2973 plasmid (Fu et al 2014) and plate on LB + 50 µg/mL Carbenicillin

#### REAGENTS

BL21-Gold (DE3) Competent Cells [230132](#) by [Agilent Technologies](#)

### Step 2.

Inoculate 25mL LB + 50 µg/mL Carbenicillin with bacteria from the fresh transformation and incubate at 37 °C overnight.

#### DURATION

18:00:00

### Step 3.

Transfer 5mL of overnight culture to 1L LB + 0.1% glucose + 50 µg/mL Carbenicillin.

### Step 4.

Grow at 25 °C to OD<sub>600</sub>=0.5.

### Step 5.

Shift culture to 18 °C for 15-25 minutes.

 **DURATION**

00:15:00

**Step 6.**

Add IPTG to 0.2 mM and incubate overnight.

 **DURATION**

18:00:00

**Step 7.**

Pellet culture and obtain wet weight.

**Step 8.**

Resuspend at 6 mL/g cells with Buffer A + protease inhibitor (Roche, #11836170001) + 1mM PMSF.

 **REAGENTS**

cComplete™, Mini, EDTA-free (Protease Inhibitor) [#11836170001](#) by [Roche](#)

**Step 9.**

Sonicate 6 x 45s (setting 3 at 30%, 1 second pulse-2 second pause) with 1 minute cooling in between.

**Step 10.**

Spin lysate 30 minutes at 16000xg and transfer supernatant to a fresh tube.

 **DURATION**

00:30:00

**Step 11.**

Equilibrate a 5mL Ni-agarose (Qiagen, #30410) with column with Buffer A (at least 25mL).

 **REAGENTS**

Ni-NTA Superflow [30410](#) by [Qiagen](#)

**Step 12.**

Batch bind clarified lysate with Ni-agarose 45 minutes at 4 °C.

 **DURATION**

00:45:00

**Step 13.**

Wash Ni-agarose column with 100mL of Buffer B.

**Step 14.**

Elute protein with Buffer C.

**Step 15.**

Determine fractions that have Cas9 protein using Bradford assay or by running a small amount on SDS-PAGE gel. Pool fractions.

Remove contaminating DNA in the prep

**Step 16.**

Equilibrate a 5mL Q Sepharose (Sigma, #Q1126) column with 1M KCl (25mL, this charges the column).

 **AMOUNT**

5 ml Additional info:

 **REAGENTS**

Q Sepharose® Fast Flow [Q1126](#) by [Sigma Aldrich](#)

Remove contaminating DNA in the prep

**Step 17.**

Equilibrate Q Sepharose column with Buffer C (25mL).

### Step 18.

Flow eluent (from step 17) over Q Sepharose column. Collect flow-through and dialyze into 1L Buffer D for 5 hours at 4 °C.

 DURATION

05:00:00

### Step 19.

Transfer into 1L Buffer D and dialyze overnight.

 DURATION

18:00:00

### Step 20.

Concentrate protein to 10 mg/mL using a 100K centrifugal filter (Milipore, UFC910024).

 REAGENTS

Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-100 membrane [UFC910024](#) by [Emd Millipore](#)

### Step 21.

Aliquot and flash-freeze in liquid nitrogen. Store aliquots at -80°C.

 NOTES

**Alexandre Paix** 04 Sep 2015

Typical yield is sufficient for 50-70 single-use aliquots (5µl aliquot, 10µg/µl Cas9).