

# Construction of PLC cells stably overexpressing Csseverin

# Mengchen Shi

## **Abstract**

Citation: Mengchen Shi Construction of PLC cells stably overexpressing Csseverin. protocols.io

dx.doi.org/10.17504/protocols.io.kcdcss6

Published: 18 Oct 2017

## **Protocol**

## Step 1.

## Step 2.

The PEZ-LV203 vector and Csseverin gene fragments were digested with EcoRI and Apa I, respectively, and subsequently ligated using T4 DNA ligase.

## Step 3.

The recombinant plasmid pEZ-LV203-Csseverin was identified by enzyme digestion and sequencing.

## Step 4.

To generate the lentivirus, the pEZ-LV203-Csseverin plasmid or PEZ-LV203 control plasmid was cotransfected into 293T cells.

### Step 5.

Supernatant containing the recombinant lentiviral particles was collected at 48 h post-transfection, filtered by a Millipore filter and subjected to ultracentrifugation.

#### Step 6.

The lentiviral particles were re-suspended in cold phosphate-buffered saline (PBS) and used to infect PLC cells.

## Step 7.

After 48 h, the cells were incubated in selection medium containing puromycin (3 mg/ml) for 7 days to select stably cells.