

# FirstStrand with RT Reaction

## New England Biolabs

### Abstract

**Citation:** New England Biolabs FirstStrand with RT Reaction. [protocols.io](https://www.protocols.io)  
[dx.doi.org/10.17504/protocols.io.cjvun5](https://doi.org/10.17504/protocols.io.cjvun5)  
**Published:** 07 Nov 2014

## Protocol

### Step 1.

RNA solution 0.5-2  $\mu\text{g}$  ( total RNA or 50-100 ng polyA-selected RNA)

### Step 2.

Primer (dT or N9) **2  $\mu\text{L}$**

### Step 3.

dNTP mix **4  $\mu\text{L}$**

### Step 4.

nuclease-free H<sub>2</sub>O to final volume of **16  $\mu\text{L}$**