

# Cas9/sgRNA ribonucleoprotein nucleofection using Lonza 4D nucleofector

Bao Thai

## Abstract

**Citation:** Bao Thai Cas9/sgRNA ribonucleoprotein nucleofection using Lonza 4D nucleofector. **protocols.io**  
dx.doi.org/10.17504/protocols.io.mpf5jn

**Published:** 17 Jan 2018

## Protocol

### Prepare cells (part 1)

#### Step 1.

Trypsinize cells and spin down at 100 x g for 5 minutes.

### Prepare cells (part 1)

#### Step 2.

Remove trypsin and resuspend cells in an appropriate amount of fresh media.

### Prepare cells (part 1)

#### Step 3.

Count cells. Record the cell concentration (cells/uL). In the meantime, put media containing cells in a 37C water bath.

### Prepare ribonucleoproteins (RNPs) mix

#### Step 4.

Add 2.5 uL of 40mM Cas9 (100 pmol) to 2.5 uL of Cas9 buffer (20 mM HEPES-KOH pH 7.5, 150 mM KCl, 10% glycerol, 1 mM TCEP-can make this ahead of time, aliquot and store at -20C).

### Prepare ribonucleoproteins (RNPs) mix

#### Step 5.

Add 3880 ng of sgRNA (120 pmol, MW32,327g/mol) to Cas9 buffer totaling 5 uL.

Prepare ribonucleoproteins (RNPs) mix

**Step 6.**

Add Cas9 to sgRNA slowly while swirling pipette tip.

Prepare ribonucleoproteins (RNPs) mix

**Step 7.**

Incubate at 37C for 10-20 minutes to let RNP form.

Prepare cells (part 2)

**Step 8.**

For each nucleofection, pipette 200k cells using a P200 or larger into a 1.5 mL tube.

Prepare cells (part 2)

**Step 9.**

Spin 100 x g for 10 minutes to pellet cells softly.

Prepare cells (part 2)

**Step 10.**

While the cells are spinning, prepare a 12-well plate containing 1 mL of media per well. Pre-warm at 37C.

Nucleofection

**Step 11.**

Prepare and label wells on nucleofection cuvettes. To avoid cells staying in nucleofection solution for a long period of time in the subsequent steps, configure Lonza 4D ahead of time using the recommended cell-type program. Use SF cell line program CM-130 for HEK293T cells.



**REAGENTS**

Lonza Nucleofector 4d [AAF-1002X](#) by [Lonza](#)

Amaza SF Cell Line 4D-Nucleofector Kit S (96 RCT) V4SC-2096 by [Lonza](#)

Nucleofection

**Step 12.**

After centrifugation, cell pellets are soft so carefully remove media from cells.

Nucleofection

**Step 13.**

Resuspend cells in 20 uL of nucleofector solution (SF cell line solution with added supplement for

HEK293T) using a P200.

#### Nucleofection

##### **Step 14.**

Add the entire 10 uL RNP mix to the 20 µL resuspension and mix using a P200.

#### Nucleofection

##### **Step 15.**

If using a repair template, add 1uL of 100uM single-stranded donor DNA (100 pmoles) and mix well.

#### Nucleofection

##### **Step 16.**

Add nucleofection mixes to the multiwell cuvette, and cap.

#### Nucleofection

##### **Step 17.**

Insert cuvette into nucleofector and zap using the configured program.

#### Nucleofection

##### **Step 18.**

Allow cells to sit in nucleofection strips for 10 minutes post-nucleofection. This is supposed to increase efficiency.

#### Nucleofection

##### **Step 19.**

Add 80uL of pre-warmed media to each well. Pipette mixture out with a P200 into your pre-warmed 12-well plate.

#### Nucleofection

##### **Step 20.**

Allow cells 24 hours to settle and recover before attempted downstream analysis. Consider including un-zapped controls to test viability.