

Column-free plasmid miniprep

Vladimir Vigdorovich

Abstract

This is a protocol for plasmid minipreps that does not involve any columns and uses commonly accessible reagents.

While it doesn't result in highly pure DNA, it can provide a sufficient quantity of plasmid for sequencing, restriction analysis, preparative digests.

Beware of the large amount of RNA that co-purifies with plasmid DNA (it acts as an efficient carrier during the nucleic acid precipitation step). This contamination will make it impossible to determine DNA concentration from standard 280-nm absorbance readings.

Credits: I first learned this protocol from the old-style protocols for bacmid purification using the Bak-to-Bac system (Invitrogen).

In that system, the large (>50 kb) bacmid molecules were deemed too fragile for the Qiagen spin columns.

Citation: Vladimir Vigdorovich Column-free plasmid miniprep. **protocols.io**

[dx.doi.org/10.17504/protocols.io.d7e9jd](https://doi.org/10.17504/protocols.io.d7e9jd)

Published: 02 Dec 2016

Guidelines

Qiagen plasmid purification buffers (from [OpenWetWare.org](https://openwetware.org))

Buffer P1

- * 50 mM Tris-HCl pH 8.0
- * 10 mM EDTA
- * 100 ug/ml RNaseA

The buffer and RNaseA can also be ordered from Qiagen separately (catalog numbers 19051 and 19101).

Note: LyseBlue reagent (1000X = 43 mg/ml thymolphthalein in ethanol)

Buffer P2

* 200 mM NaOH

* 1% SDS

Buffer P3 (not for spin columns, but for Qiatips, midi, maxi, giga kits; N3 contains the equilibration buffer, while P3 does not)

* 3.0 M potassium acetate pH 5.5

Materials

- ✓ Buffer P1 by Contributed by users
- ✓ Buffer P2 by Contributed by users
- ✓ Buffer P3 19053 by Contributed by users
- ✓ Isopropanol by Contributed by users
- ✓ 1.5 mL Eppendorf tubes by Contributed by users
- ✓ nuclease-free water by Contributed by users

Protocol

Culture growth

Step 1.

Grow an overnight 3-mL *E. coli* culture in LB with appropriate antibiotics. Preferably, in round-bottom snap-cap culture tubes suitable for centrifugation (in the next step).

AMOUNT

3 µl Additional info:

Culture harvesting

Step 2.

Spin down at 4,000 rpm in a table-top centrifuge for 15 min. Discard supernatant and vortex pellets to resuspend in residual culture medium.

DURATION

00:15:00

Resuspension

Step 3.

Add the roughly resuspended cells to 200 μ L Buffer P1. Mix by pipetting or vortexing.

 AMOUNT

200 μ L Additional info:

Lysis

Step 4.

Add 200 μ L Buffer P2.

Mix gently by inverting the tube several times.

 AMOUNT

200 μ L Additional info:

Neutralization

Step 5.

Add 350 μ L Buffer P3.

Mix vigorously by inverting/shaking the tube several times. Do not vortex.

Be sure to thoroughly mix the viscous solution from previous step and allow the flocculant SDS precipitate to form.

 AMOUNT

350 μ L Additional info:

Lysate clarification

Step 6.

Centrifuge for 5 min in a microcentrifuge ($>10,000$ rpm) to remove SDS and cell debris.

 DURATION

00:05:00

Nucleic acid precipitation

Step 7.

Transfer the clarified lysate (750 μ L) to a clean tube containing 750 μ L pure isopropanol.

Mix by inversion.

 AMOUNT

750 μ L Additional info:

Nucleic acid isolation

Step 8.

Centrifuge for 5 min in a microcentrifuge ($>10,000$ rpm) to pellet nucleic acids.

 DURATION

00:05:00

Nucleic acid pellet cleanup

Step 9.

Wash the nucleic acid pellet with 500 µL of 70% ethanol.

Gently dislodge the nucleic acid pellet from the side of the tube by shaking.

AMOUNT

500 µl Additional info:

Nucleic acid pellet cleanup: salt removal

Step 10.

Centrifuge for 5 min in a microcentrifuge (>10,000 rpm) to pellet nucleic acids.

DURATION

00:05:00

Nucleic acid pellet cleanup: ethanol removal

Step 11.

Gently pipette off the ethanol and air-dry the nucleic acid pellet.

If you have a SpeedVac centrifuge, use it for 2 min.

Don't overdry the pellets, as that will cause problems with resuspension.

Nucleic acid pellet resuspension

Step 12.

Dissolve the pellet in 50 µL nuclease-free water.