

Protocol for transfection of H4, (ATCC®: HTB-148) Cells by FuGENE HD in 96 well plates.

Misha Gurevich, V. Katerov

Abstract

This is a protocol for transfection of H4 (Human Brain Neuroglioma) in either 10% or 100% FBS using FuGENE HD.

Citation: Misha Gurevich, V. Katerov Protocol for transfection of H4, (ATCC®: HTB-148) Cells by FuGENE HD in 96 well plates.. **protocols.io**

dx.doi.org/10.17504/protocols.io.e8mbhu6

Published: 23 Jun 2016

Protocol

Step 1.

Cell plating

H4 (neuroglioma) cells were seeded from 90-100% confluent culture the day before transfection with the density 7,000 cells/well in 100µl complete growth medium (DMEM + 10% Fetal Bovine Serum).

For transfection in 100% Fetal Bovine Serum the complete growth medium was replaced with Fetal Bovine Serum two hours before transfection.

Step 2.

Complex preparation (per 20 wells)

Tissue culture 96-round bottom well plates were used for complex preparation:

Step 3.

Prepare 0.02µg/µl pCMVβ DNA solution in OptiMEM® or sterile deionized water.

Step 4.

Add 6 µl of reagent to 100 µl of DNA solution.

Step 5.

Mix carefully by pipetting (15 times).

Step 6.

Incubate 10 min at room temperature.

Step 7.

Add 5 µl of complex per well to the cells, and mix thoroughly.

Step 8.**Incubation**

Incubate transfected cells in CO₂ incubator for 48 hours.

Step 9.**Detection of β-gal expression**

1. Remove the medium from the well and wash the cells once with 100µl per well PBS.

Step 10.

Fix the cells in the well with 50µl solution of 4% formaldehyde in PBS for 5min at room temperature.

Step 11.

Wash each well twice with 100µl PBS.

Step 12.

Add 50µl per well of substrate/stain solution and incubate the plate overnight at 37°C.

Step 13.

Observe the cells under microscope and evaluate the proportion of blue (β-gal-positive) cells.