# **8 Marine Purity Broth**

# **Bonnie Poulos**

# **Abstract**

Purpose: To test axenic cultures for purity.

Protocol described in S. Bertillson, O. Berglund, D.M. Karl, S.W. Chisholm (2003). Elemental composition of marine Prochlorococcus and Synechococcus: Implications for the ecological stoichiometry of the sea. Limnol Oceanogr 48(5):1721-1731.

Citation: Bonnie Poulos Marine Purity Broth. protocols.io

dx.doi.org/10.17504/protocols.io.ekrbcv6

Published: 23 Feb 2016

# **Guidelines**

#### **Formulation:**

17 g BD Difco AC Broth 20 g NaCl Qs to 1 L with distilled H<sub>2</sub>O

## After broth has cooled:

2 M MgSO<sub>4</sub> (2 M MgSO<sub>4</sub>·7H<sub>2</sub>O = 49.3 g/100 ml H<sub>2</sub>O, dissolve and 0.2  $\mu$ m filter) 1 M CaCl<sub>2</sub> (1 M CaCl<sub>2</sub>·2H<sub>2</sub>O = 14.7 g/100 ml H<sub>2</sub>O, dissolve and 0.2  $\mu$ m filter)

# **Protocol**

### Step 1.

Autoclave BD Difco AC Broth and NaCl Os to 1 L with distilled water for 20 min.

**O DURATION** 

00:20:00

#### **P** NOTES

Bonnie Poulos 23 Feb 2016

See formulation in guidelines.

#### Step 2.

After broth has cooled, add 16 ml 2 M MgSO<sub>4</sub> and 10 ml 1 M CaCl<sub>2</sub>

#### NOTES

Bonnie Poulos 23 Feb 2016

(2 M MgSO<sub>4</sub>·7H<sub>2</sub>O = 49.3 g/100 ml H<sub>2</sub>O, dissolve and 0.2  $\mu$ m filter)

Bonnie Poulos 23 Feb 2016

(1 M CaCl<sub>2</sub>·2H<sub>2</sub>O = 14.7 g/100 ml H<sub>2</sub>O, dissolve and 0.2  $\mu$ m filter)

#### Step 3.

Dispense 2.0 or 4.5 ml per sterile tube

# Step 4.

Inoculate sterile MPB with 1:5 (vol:vol) dilution of culture

### Step 5.

Incubate at 22°C

# NOTES

# Bonnie Poulos 23 Feb 2016

Heterotrophs normally grow in less than 2 days. No growth in the broth after 1-2 weeks indicates culture is axenic.