

# BiTE® Xenograft Protocol

Sandra L. Ross, Marika Sherman, Patricia L. McElroy, Julie A. Lofgren, Gordon Moody, Patrick A. Baeuerle, Angela Coxon, Tara Arvedson

## Abstract

This protocol describes methods for a mouse xenograft study designed to measure in vivo BiTE® efficacy against antigen-expressing human tumor cells. In a modified version, the protocol is used to measure BiTE®-mediated bystander killing. Mixtures of human EGFR-positive (unlabeled) and EGFR-negative (luciferase-labeled) cells are implanted in nude mice along with human T cells. Mice are treated daily with BiTE®, beginning one day after tumor implant for 20 days. Growth of EGFR-negative cells is measured by imaging (luminescence) and by calipers (tumor volume).

**Citation:** Sandra L. Ross, Marika Sherman, Patricia L. McElroy, Julie A. Lofgren, Gordon Moody, Patrick A. Baeuerle, Angela Coxon, Tara Arvedson BiTE® Xenograft Protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.h5eb83e

**Published:** 22 May 2017

## Protocol

### Overview

#### Step 1.

- Female athymic nude mice approximately 7 weeks old (Charles River Laboratories) are used for the study
- Human tumor cell lines are cultured in McCoy's medium with 10% heat-inactivated fetal bovine serum (FBS) (HCT116) or RPMI with 10% FBS (SW620-LUC) at 37° C with 5% CO<sub>2</sub>.
- Human pan T cells (AllCells) are thawed according to supplier instructions and cultured in RPMI with 10% FBS at 37° with 5% CO<sub>2</sub>. T cells are activated by culturing with beads coated with anti-CD2, -CD28 and CD3 antibodies following the manufacturer's instructions (Miltenyi Biotec T cell activation/expansion kit). Three days after establishment of the culture, hIL-2 (Miltenyi Biotec) is added to the medium at a concentration of 0.1 µg/mL.
- Tumor cells mixed with activated T cells at a 1:1 ratio are implanted subcutaneously on the left flank of the mice, in a 1:1 mix of serum-free RPMI and Matrigel basement membrane matrix (BD Biosciences).
- Control and treatment BiTEs® are diluted to the proper concentration in PBS and administered intraperitoneally.
- Imaging of bioluminescence from luciferase-labeled cells is conducted using the IVIS Spectrum in vivo imaging system (PerkinElmer); mice are administered 150 mg/kg D- Luciferin (Sigma, diluted in PBS) intraperitoneally 15 minutes before image capture.
- Tumor volumes are obtained using handheld digital calipers to measure tumor length, width and height to calculate cubic millimeters; mice are weighed using a digital balance.
- Tumors are collected at the end of the study and placed in neutral buffered formalin (NBF).

### Prepare T cells

## Step 2.

1. Prepare activating beads according to manufacturer's instructions; briefly, prepare mixture of activating antibodies, add unlabeled beads, mix gently at 4° for 2 hours, store at 4°.
2. Thaw T cells according to supplier instructions and place in tissue culture (TC) flask, at a concentration of approximately  $2 \times 10^6$  cells/mL.
3. Add activating beads at a bead-to-cell ratio of approximately 1:4. Wash beads one time to remove excess antibody before adding to T cells.
4. After 3 days, add medium with IL-2. Keep cell density at approximately  $2 \times 10^6$  cells/mL, adding fresh media with IL-2 and expanding into additional flasks every 2 to 3 days. Viability should remain  $\geq 90\%$ . Expand to approximately 50% more than number needed for implantation; extra is needed due to loss during de-beaded prior to implantation.

## Prepare tumor (target) cells

### Step 3.

1. Thaw tumor cells and place in proper medium in TC flasks and incubate.
2. Expand to cell numbers necessary for implant (plan for approximately 30% extra cells for each cell line).

## Implant admixed tumor (target) cells and T cells

### Step 4.

Day 0:

1. Thaw a vial of Matrigel; keep cold.
2. Count T cells to determine viability and cells/mL. Viability should be  $\geq 90\%$ . Determine volume of T cells needed (with extra) based on cell count and remove beads by using MACSiMAG Separator magnetic column (Miltenyi Biotec) or similar. Count de-beaded T cells and keep on ice.
3. Harvest appropriate number of tumor cells using 0.05% trypsin with EDTA (1X) (Gibco), adding medium with FBS after cells have lifted off the flask; pool cells, determine cell count and viability (should be  $\geq 90\%$ ); keep on ice.
4. Mix T cells and tumor cells at a 1:1 ratio. (Example: add 250,000 T cells to 250,000 tumor cells.) For groups containing two tumor cell lines, combine equal numbers of each tumor cell line (e.g. 250,000 HCT116 and 250,000 SW620-LUC) and add T cells equal to total number of tumor cells (500,000).
5. Wash T cell and tumor cell mixture in serum-free RPMI (centrifuge at 300 G, 5 minutes, 4°). Resuspend pellet in a 1:1 mixture of serum-free RPMI and Matrigel, in a volume that will give each mouse the desired number of cells in 0.1 mL. Keep cell solution on ice.
6. Inject 0.1 mL of T cell/tumor cell mixture subcutaneously on the flank of each mouse, using a 0.5 mL syringe with a 28 gauge needle. Fill only one syringe at a time and keep remaining cell solution on ice so that the Matrigel does not thicken.

## Baseline bioluminescence and initiation of BiTE® treatment

### Step 5.

Day 1:

1. Image all mice with IVIS bioluminescence imaging system to obtain baseline image. Treat mice with 150 mg/kg IP of D-Luciferin 15 minutes prior to imaging. Image tumors with mice in lateral

recumbency with tumor side up.

2. Dilute BiTEs® with PBS to proper concentration.
3. Treat mice with control or active BiTE®: administer IP.

#### Ongoing BiTE® dosing and tumor measurement

##### **Step 6.**

Day 2 and continuing until end of study:

- Treat mice daily with BiTEs®.
- Image 2 times per week, and measure tumors with calipers.
- Weigh mice 2 times per week once tumors are palpable.
- At end of study, euthanize mice according to IACUC guidelines and harvest tumors; place in NBF.