

# First Strand Synthesis with Reverse Transcriptase

## New England Biolabs

### Abstract

This protocol is for First Strand Synthesis with Reverse Transcriptase

**Citation:** New England Biolabs First Strand Synthesis with Reverse Transcriptase. **protocols.io**

dx.doi.org/10.17504/protocols.io.cjuunv

**Published:** 27 Jan 2015

## Protocol

### Step 1.

In a sterile microfuge tube add the following:

 [PROTOCOL](#)

#### . [FirstStrand with RT Reaction](#)

CONTACT: [New England Biolabs](#)

#### Step 1.1.

RNA solution 0.5-2 µg ( total RNA or 50-100 ng polyA-selected RNA)

#### Step 1.2.

Primer (dT or N9) **2 µL**

#### Step 1.3.

dNTP mix **4 µL**

#### Step 1.4.

nuclease-free H<sub>2</sub>O to final volume of **16 µL**

### Step 2.

Heat for 3-5 minutes at 65-80°C.

 [DURATION](#)

00:05:00

### Step 3.

Spin briefly and place promptly on ice.

### Step 4.

Add 10X RT Buffer

 [AMOUNT](#)

2 µl Additional info:

 [PROTOCOL](#)

#### . [10X RT Buffer](#)

CONTACT: [New England Biolabs](#)

**Step 4.1.**

500 mM Tris-HCl (pH 8.3 @ 25°C)

**Step 4.2.**

750 mM KCl

**Step 4.3.**

30 mM MgCl<sub>2</sub>

**Step 4.4.**

100 mM DTT

**Step 5.**

Add RNase inhibitor **1 µL**

**Step 6.**

Add M-MuLV Reverse Transcriptase **1 µL**

**Step 7.**

Incubate at 42°C for one hour.

 **DURATION**

01:00:00

**Step 8.**

Inactivate enzyme at 90°C for 10 minutes.

 **DURATION**

00:10:00

**Step 9.**

Store products at -20°C or proceed to next step(s).