

# Viral DNA Miniprep Procedure

David Dunigan and Irina Agarkova

## Abstract

**Citation:** David Dunigan and Irina Agarkova Viral DNA Miniprep Procedure. **protocols.io**

dx.doi.org/10.17504/protocols.io.eq6bdze

**Published:** 13 Jun 2016

## Guidelines

### Materials:

- 1) 60-65°C heat block or water bath
- 2) Microfuge
- 3) 1.5 and 2.0 mL microfuge tubes (screw-cap)
- 4) 50 mM Tris-HCl, pH 7.5, 10 mM MgCl<sub>2</sub>
- 5) Triton X-100
- 6) DNase I, 2.0 mg/mL in 50 mM Tris-HCl, pH 8.0. Store in 110 µL aliquots at -20°C DO NOT REFREEZE UNUSED MATERIAL, DISCARD.
- 7) Proteinase K, 2.0 mg/mL in 50 mM Tris-HCl, pH 8.0. Autodigest for 60 min at 37°C before use. Store in 1.0 mL aliquots at -20°C. Can be refrozen unless the material has fallen out of solution.
- 8) 10% Na sarcosyl
- 9) CHCl<sub>3</sub>:Isoamyl alcohol (24:1)
- 10) 500 mM EDTA, pH 8.0
- 11) 3 M NaOAc
- 12) 100% EtOH
- 13) Buffer-saturated phenol
  - Preparation: Thaw 100 gm bottles of phenol at 60-65°C. Add a stir bar, 100 mL of 500 mM Tris-HCl, pH 8.0 and 0.1 gm 8-hydroxyquinoline. Stir and allow the phases to separate at 4°C overnight. Aspirate off the upper aqueous layer and add 75 ml of 100 mM Tris-HCl, pH 8.0, 0.2% 2-mercaptoethanol (2-ME) and stir. Allow the phases to separate at 4°C for several hours to overnight and remove the upper aqueous layer. Repeat the 75 mL addition 2X, leaving the final

phase on the phenol. Store at 4°C.

14) 10 mM Tris-HCl, pH 8.0, 1 mM EDTA (1X TE)

## Protocol

### Step 1.

Infect 60 mL of chlorella with 200 µL of viral single plaque isolates.

### Step 2.

Incubate the samples at 25°C for 24-72 hours, with continuous light and shaking.

 DURATION

12:00:00

### Step 3.

Centrifuge 30 mL of the lysates in the Sorvall SS34 rotor at 5,000 rpm (3,000 rcf), 5 min, 4°C.

 DURATION

00:05:00

### Step 4.

Save the supernatants. Save the unused portion of the lysates.

### Step 5.

Add 10% NP-40 (or Triton X-100) to the lysate supernatants to a final concentration of 1%.

### Step 6.

Centrifuge the material in Beckman Ti50.2 rotors at 15,000 rpm (27,000 rcfmax), 75 min, 4°C.

 DURATION

01:15:00

### Step 7.

Discard the supernatants.

### Step 8.

Resuspend the virus pellets with 1.0 mL of 50 mM Tris-HCl, pH 7.5, 10 mM MgCl<sub>2</sub>.

### Step 9.

Transfer 350 µL of the resuspended virus to 1.5 mL screw-cap microfuge tubes and adjust the final volume to 500 µL with 50 mM Tris-HCl, pH 7.5, 10 mM MgCl<sub>2</sub>.

### Step 10.

Add 8.8 µL of DNase I and mix.

**Step 11.**

Incubate at room temperature for 60 min.

 DURATION

01:00:00

**Step 12.**

Add 6.0 µL of 500 mM EDTA, pH 8.0 to the samples and mix.

**Step 13.**

Add 56.6 µL of proteinase K and 29.0 µL of 10% Na sarcosyl and mix.

**Step 14.**

Incubate the samples at 60-65°C for 60 min.

 DURATION

01:00:00

**Step 15.**

Add 300 µL of buffer-saturated phenol and 300 µL of CHCl<sub>3</sub>:Isoamyl alcohol (24:1) to the tubes.

**Step 16.**

Mix by inversion.

**Step 17.**

Centrifuge in the microfuge at maximum speed for 5 min at 4°C.

 DURATION

00:05:00

**Step 18.**

Remove the upper aqueous layers to clean tubes.

**Step 19.**

Add 600 µL of CHCl<sub>3</sub>:Isoamyl alcohol (24:1) to the tubes.

**Step 20.**

Mix by inversion and centrifuge for 5 min at 4°C in the microfuge.

 DURATION

00:05:00

**Step 21.**

Remove the upper aqueous layers to clean tubes and repeat the CHCl<sub>3</sub>:Isoamyl alcohol extraction 1X.

**Step 22.**

Place the last extraction into 2.0 mL microfuge tubes.

**Step 23.**

Add 66 µL of 3 M NaOAc to each tube.

**Step 24.**

Precipitate the DNAs with 2X volumes (approximately 1350 µL) of 100% EtOH.

**Step 25.**

Mix well and hold at -20°C overnight.

 DURATION

18:00:00

**Step 26.**

Centrifuge the tubes in the microfuge for 10-15 min at 4°C to pellet the DNAs.

 DURATION

00:15:00

**Step 27.**

Discard the supernatants.

**Step 28.**

Wash the DNA pellets 1X with 1000  $\mu$ L of 70% EtOH in the microfuge for 5 min at 4°C.

 **DURATION**

00:05:00

**Step 29.**

Dry the pellets briefly (10-15 min) in the vacuum desiccator or the speed vac (5 min) to remove the EtOH.

 **DURATION**

00:15:00

**Step 30.**

Resuspend the DNAs with approximately 60  $\mu$ L of 1X TE buffer. If the DNA doesn't go into solution overnight, centrifuge in the microfuge for 15 min at 4°C and remove the supernatants to clean tubes.

**Step 31.**

Discard the pellets.

**Step 32.**

Store the DNAs at 4°C.