

May 09, 2019

Working

## U Mass - Alanine Transferase 👄

Jason Kim<sup>1</sup>

<sup>1</sup>University of Massachusetts

dx.doi.org/10.17504/protocols.io.xpmfmk6

Mouse Metabolic Phenotyping Centers

 $Tech. \ support\ email: \textbf{info@mmpc.org}$ 



Lili Liang 69

ABSTRACT

### Summary:

This experiment involves a spectrophotometric measurement using Roche Cobas Clinical Chemistry Analyzer. Serum alanine transferase levels reflect liver function.

**EXTERNAL LINK** 

https://mmpc.org/shared/document.aspx?id=180&docType=Protocol

#### MATERIALS

NAME ~	CATALOG #	VENDOR ~	CAS NUMBER $\vee$ RRID $\vee$
Pyridoxal Phosphate	04774221190	Roche	
Calibrator f.a.s.	10759350 360	Roche	
Precinorm U Plus	12149435 160	Roche	
Precipath U Plus	12149443 160	Roche	
NaCl Diluent 9%	04774230 190	Roche	
Cleaner	04774248 190	Roche	
Micro Sample cups	11406680 001	Roche	
Chimneys	11930630 001	Roche	
NERL High Quality Water	9805	Fisher Scientific	
Alanine aminotransferase acc. IFCC	04718569190	Roche	

MATERIALS TEXT

# Note:

Roche, RRID:SCR\_001326

Fisher Scientific, RRID:SCR\_008452

BEFORE STARTING

#### Notes:

- $\,\,\sqrt{\,}\,$  Try to use freshly prepared serum and plasma samples for this assay.
- √ No dilution or treatment of the sample is required, but plasma samples should be centrifuged to remove any fibrin/fibrinogen
- √ Samples should be stored at 2-8°C for 24 hours prior to analysis. For longer periods, store samples at -70°C, and avoid repeated

freeze/thaw cycle
-------------------

√ A 50 µldead volume is required in addition to sample volume for multi-protein analysis (typically 1-5 µl).

- 1 Perform daily quality control assessment of instrumentation before analysis.
- 2 Load each sample into a specialized micro-sample cup for the clinical chemistry analyzer.
- 3 Select ALT test on display and run the analysis.
- 4 Collect and analyze the data.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited