



Apr 17, 2019

In devel.

## ActA purification protocol

Sharifah Albraiki<sup>1</sup><sup>1</sup>Wichita State University[dx.doi.org/10.17504/protocols.io.nkadcse](https://doi.org/10.17504/protocols.io.nkadcse)

Beck Lab

Sharifah Albraiki  
Wichita State University

### Buffers

- Lysis Buffer:**  
20mM MOPS  
100mM KCl  
5mM Imidazole pH 7  
  
**Wash Buffer:**  
20mM MOPS  
100mM KCl  
20mM Imidazole pH 7  
  
**Elution Buffer:**  
20mM MOPS  
100mM KCl  
250mM Imidazole pH 7  
  
**Storage Buffer:**  
20mM HEPES  
100 NaCl pH 7.5

### Day 1

- Spread the bacteria stock from -80 (tray D – first box tube name “Listeria ActA-His DPL 1545”) in BHI plate (no antibiotic), and incubate @ 37°C for overnight.
- Make Brain Heart Infusion media (BHI)

### Day 2

- Start 25ml culture + 7.14ul chloramphenicol @37°C for overnight

### Day 3

- Transfer the 25ml starting culture into (475ml BHI+135.7ul chloramphenicol) and grow @37°C for 6-8 hours
- Spin down cells in 350ml bottles at (4500xg=6000rpm) for 20-30 minutes. (save supernatant)

- 7 ■ Precipitate the protein from the supernatant by adding slowly the ammonium sulfate (measured and ground to powder) to 40% (ActA-His) or 60% (ActA-N-His) saturation. (stir in cold room) for 1 hour.
- 8 ■ Centrifuge and pellet the precipitated protein at 4500xg for 20-30 minutes.
- 9 ■ Re-suspend the precipitated protein in the lysis buffer

#### Day 3-5

- 10 ■ Dialyze the protein in the storage buffer for 3-4 times every six hours to remove the salt before using the Nickel column

#### Day 6

- 11 ■ Purify with Ni-NTA (wash 2 times and elute the protein with the elution buffer). Then, concentrate the elution.

#### Day 7

- 12 ■ Purify the eluted concentrated proteins by Superose-6 (Pharmacia) gel-filtration chromatography.(optional)



When I did the gel-filtration step I lost all my protein.

- 13 ■ freeze in N<sub>2</sub>(l) and stored at 80°C (optional)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited