Measurement of cellular oxygen consumption and extracellular acidification 👄

PLOS One

Apr 16, 2019

Working

Kiichi Hirota¹, Yoshiyuki Matsuo¹

¹Kansai Medical University

dx.doi.org/10.17504/protocols.io.v92e98e



Yoshiyuki Matsuo

EXTERNAL LINK

https://doi.org/10.1371/journal.pone.0215072

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Sumi C, Matsuo Y, Kusunoki M, Shoji T, Uba T, Iwai T, Bono H, Hirota K (2019) Cancerous phenotypes associated with hypoxiainducible factors are not influenced by the volatile anesthetic isoflurane in renal cell carcinoma. PLoS ONE 14(4): e0215072. doi: 10.1371/journal.pone.0215072

BEFORE STARTING

The cellular oxygen consumption rate (OCR) and extracellular acidification rate (ECAR) were measured using an XFp Extracellular Flux Analyzer™ (Agilent Technologies, Santa Clara, CA).

- Cells were seeded on to the XFp Cell Culture microplate and incubated overnight.
- The sensor cartridge of the XFp Analyzer was hydrated at 37 °C in a non-CO2 incubator one day before the experiment.
- For the OCR assay, injection port A on the sensor cartridge was loaded with oligomycin (a complex V inhibitor, final concentration 1 µM), 3 carbonyl cyanide-p-trifluoromethoxyphenylhydrazone (FCCP; an uncoupling agent, final concentration 2 µM) was loaded to port B, and rotenone/antimycin A (inhibitors of complexes I and III, final concentration 0.5 μM each) was loaded to port C.
- During sensor calibration, cells were incubated at 37 °C in the non-CO2 incubator in 180 µl assay medium (XF Base Medium: 25 mM glucose, 1 mM pyruvate, and 2 mM L-glutamine, pH 7.4).
- The plate was immediately placed on the calibrated XFp Analyzer for the assay.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited