

# PhosphoThreonine Western Blotting

**Christopher Bartley**

## Abstract

This protocol is for pSerine and pThreonine western blotting (optimized for detecting phospho FMRP protein after IP)

**Citation:** Christopher Bartley PhosphoThreonine Western Blotting. **protocols.io**

dx.doi.org/10.17504/protocols.io.c45yy5

**Published:** 01 Nov 2016

## Protocol

### Step 1.

Run 40% of IP eluate on 10% gel Tris-Glycine



#### REAGENTS

Anti-Phosphoserine Antibody [AB1603](#) by [Emd Millipore](#)

### Step 2.

Transfer protein to PVDF membrane for 1hr at RT with ice pack (100V)

### Step 3.

Cut membrane just below 180kD (top MW marker on benchmark prestained ladder)



#### REAGENTS

Bovine Serum Albumin (IgG-Free, Protease-Free) [001-000-161](#) by [Jackson ImmunoResearch](#)

### Step 4.

Cut Membrane at 82kD (the blue band above the pink band using benchmark prestained ladder).

**\*\*Note,** the GST-tagged protein is 115 kD. Cutting the membrane in the places will help reduce non-specific binding.

Also, running a ladder in the middle of the gel will help you cut straight across the membrane in the event that the gel is transferred at an angle relative to the membrane.



#### REAGENTS

Bovine Serum Albumin (IgG-Free, Protease-Free) [001-000-161](#) by [Jackson ImmunoResearch](#)

### Step 5.

Block membrane with 5% BSA (IgG and Protease-Free)/TBST for 2hrs at RT



#### REAGENTS

Bovine Serum Albumin (IgG-Free, Protease-Free) [001-000-161](#) by [Jackson ImmunoResearch](#)

### Step 6.

Probe membrane with phospho-threonine-HRP antibody in IgG-Free/Protease-Free BSA/TBST (1:1000). Probe overnight at 4C in cold room on tilting tray.

#### REAGENTS

Phospho-Threonine Antibody (P-Thr-Polyclonal) (HRP Conjugate) #6949 [6949](#) by [Cell Signaling Technology](#)

#### **Step 7.**

Rinse membrane 8x over the course of 2hrs with 1x TBST

#### **Step 8.**

Activate membrane with pico ECL

**\*\*Note,** the anti-pThr antibody is HRP conjugated so you don't need to use a secondary.

#### REAGENTS

✓ SuperSignal™ West Pico Chemiluminescent Substrate (Pico) [34080](#) by Contributed by users