

Transformation of *Skeletonema marinoi* using Multipulse Electroporation

Oskar N. Johansson., Adrian K. Clarke

Abstract

*The following transformation protocol is designed for the insertion of linear DNA constructs into the nuclear genome of *Skeletonema marinoi* by non-homologous recombination.*

Citation: Oskar N. Johansson., Adrian K. Clarke Transformation of *Skeletonema marinoi* using Multipulse Electroporation. **protocols.io**

dx.doi.org/10.17504/protocols.io.g6nbzde

Published: 10 Feb 2017

Protocol

Cell preparation

Step 1.

- a) Grow two 400 mL cultures in Artificial Seawater with f/2+Si supplements (Growth media) for ca. one week until dense. The standard growth conditions used in our lab are 16°C, 50-70 $\mu\text{mol. photons m}^{-2} \text{ s}^{-1}$, 16 h photoperiod.
- b) Reduce the volume of the sedimented culture in each flask to 50 mL by suction. Resuspend cells by vigorous agitation to minimize chain length and then transfer them to a 50 mL centrifuge tube.
- c) Pellet cells by centrifugation (1200 x g, 5 min, 4°C, swing-out rotor). Decant supernatant and resuspend the pellet in each tube in 10 mL of ice-cold 0.3M sorbitol.
- d) Pellet cells once more by centrifugation (1200 x g, 5 min, 4°C, swing-out rotor) and resuspend in 2.5 mL of ice-cold 0.3M sorbitol. Pool both resuspensions in a single tube.
- e) Store the cell suspension on ice until needed.

■ ANNOTATIONS

Jon Udell 17 Feb 2017

Test

Electroporation Equipment

Step 2.

This protocol was tested using the Gene Pulser Xcell electroporator (BioRad, USA) with the recommended 2 mm cuvettes (BioRad, USA).

Combining the cell suspension with desired DNA

Step 3.

Add 3-5 µg of the linear DNA construct to a 100 µL cell suspension in a 2 mm electroporation cuvette and leave for 3-5 min at 4°C.

NB. DNA construct should include a selectable marker such as antibiotic resistance (e.g., zeocin/bleomycin)

Electroporation

Step 4.

Ensure cells are resuspended in the cuvette prior to electroporation. Electroporation is performed by the following steps:

- a) Poring pulses (300 V, 6 pulses, 1 sec pulse Interval, 5 ms pulse length)
- b) Transfer pulses (40 pulses [10x4], 50 ms pulse length, 0.1 s pulse interval)
- c) Reverse direction of the cuvette.
- d) Repeat transfer pulses (40 pulses [10x4], 50 ms pulse length, 0.1 s pulse interval)

Cell recovery

Step 5.

Transfer the electroporated cells from the cuvette using a Pasteur pipette to 30 mL liquid growth media (without selection) and leave for 48 h under standard growth conditions

Apply Selection

Step 6.

Add selection to the media and continue growth for another 48 h

Transfer to Solid media

Step 7.

- a) Pellet cells in 50 mL centrifuge tubes (1200 x g, 5 min, room temperature) and resuspend in 1 mL growth medium with selection.
- b) Spread 75 (for 7 cm diameter plates) or 250 µL (for 15 cm diameter plates) of the resuspension on growth media plates (f/2 + Si + selection, 0.9% Agar).

Left under standard growth conditions, colonies should appear after 2-3 weeks.