

# IMPACT purification of native proteins Version 2

**Brian Baltzar**

## Abstract

**Citation:** Brian Baltzar IMPACT purification of native proteins. **protocols.io**

[dx.doi.org/10.17504/protocols.io.gj2buqe](https://dx.doi.org/10.17504/protocols.io.gj2buqe)

**Published:** 28 Nov 2016

## Guidelines

### IMPACT Column Buffer (1L)

Reagent	Volume
Tris-HCl 7-8 pH 1M	20 ml
NaCl 5M	100 ml
EDTA 0.5M	2 ml
Tween-20	2 ml
H <sub>2</sub> O	-> 1000ml

### IMPACT Cleavage Buffer (100ml)

Reagent	Volume
Tris-HCl 7-8 pH 1M	2 ml
NaCl 5M	10 ml
EDTA 0.5M	0.2 ml
DTT	0.2 ml
H <sub>2</sub> O	-> 100ml

## Before start

Make sure that the desired protein has the correct restriction sites and is in frame with the intein in the purification vector.

1. If that is not the case natively, design primers using the sequences in the following table.

RESTRICTION SITE	SEQUENCE (RESTRICTION SITE UNDERLINED)	CLONING VECTOR
NdeI	5'- GGT GGT <u>CAT ATG</u> NNN NNN... -3' (forward primer)	pTXB1
SapI <sup>1</sup>	5'- GGT GGT <u>TGC TCT TCC</u> GCA NNN NNN...-3' (reverse primer)	pTXB1
SapI <sup>2</sup>	5'- GGT GGT <u>TGC TCT TCC</u> AAC NNN NNN... -3' (forward primer)	pTYB21
PstI <sup>3</sup>	5'- GGT GGT <u>CTG CAG</u> TCA NNN NNN... -3' (reverse primer)	pTYB21

<sup>1</sup> SapI digestion creates a 3-nt overhang (GCA) for ligation with the SapI-digested pTXB1 vector (containing a TGC overhang), resulting in an in-frame fusion to the N-terminus of an intein. The SapI site can be used to add one or more extra amino acid residue(s) to the target protein by including an appropriate sequence (e.g. add ACC in the reverse primer corresponding to a GGT codon for a glycine residue). The SapI site is not regenerated after cloning.

<sup>2</sup> SapI digestion creates a 3-nt overhang (AAC) compatible with the SapI digested pTYB21 (containing a GTT overhang). The SapI site is not regenerated after cloning.

<sup>3</sup> A stop codon should be included in the reverse primer when constructing a N-terminal fusion.

2.

3. Then perform a Phusion PCR with the designed primers

## Materials

🐛 Chitin Resin - 100 ml [S6651L](#) by [New England Biolabs](#)

EDTA (0.5 M), pH 8.0 [AM9260G](#) by [Life Technologies](#)

✓ Sodium Chloride [PubChem CID: 5234](#) by Contributed by users

Tween 20 [170-6606-MSDS](#) by [Bio-rad Laboratories](#)

Econo-Pac Columns [7321010](#) by [Bio-rad Laboratories](#)

## Protocol

Insert gene into IMPACT vector

### Step 1.

Digest the pPCR product with NdeI & SapI in minimum 2 hours.

### Step 2.

Purify the digest on columns

### Step 3.

Ligate the purified digest products using the preferred ligation protocol

**Step 4.**

Make a TSB transformation into *E. Coli*:K-12 use your preferred transformation protocol

**Step 5.**

Inspect the ligation with a colonyPCR

**Step 6.**

Make ON of the verified transformation

**Step 7.**

Purify the plasmid from the ON from step 1.7 using Miniprep

**Protein expression****Step 8.**

Transform the minipreped plasmid from step 7 into *E. Coli*:ER2566

**Step 9.**

Inoculate 1L LB+Amp with a fresh colony from the transformation above.

**Step 10.**

Incubate at 37°C till it reaches a OD<sub>600</sub> at 0.5-0.8

**Step 11.**

Add IPTG for a end concentration of 0.4 mM

**Step 12.**

Incubate at 16-20°C overnight.

**Step 13.**

Centrifuge the 1L ON at 5000G at 4°C - remove the supernatant.



00:15:00

#### Step 14.

Resuspend the cellpellet in 100 ml cold column buffer - (See recipe in description)

 AMOUNT

100 ml Additional info:

#### Step 15.

Centrifuge at 8000G at 2°C - remove the supernatant.

 DURATION

00:20:00

#### Step 16.

Resuspend the cellpellet in a small volume of cold column buffer (2-4 mL)

#### Step 17.

Lyse the cells using french press

#### Step 18.

Centrifuge the lysant at 15000G at 4°C

 DURATION

00:30:00

Binding to Chitin beads

#### Step 19.

Take 10 ml Chitin beads put into a 50ml falcon tube

 AMOUNT

10 ml Additional info:

#### Step 20.

Centrifuge at 200G - remove the supernatant

 DURATION

00:05:00

#### Step 21.

Wash with 50 ml Column Buffer

 AMOUNT

50 ml Additional info:

**Step 22.**

Centrifuge at 200G - remove the supernatant

 DURATION

00:05:00

**Step 23.**

Load the supernatant from step 18

**Step 24.**

Attach on a rotor with slow spin

 DURATION

02:00:00

**Step 25.**

Add to a column, open valve just enough that it slowly drips out.

**Step 26.**

Add column Buffer.

 AMOUNT

200 ml Additional info:

**Step 27.**

When the column buffer has run through, close the valve and add Cleavage Buffer (just enough that it covers the beads + 1ml.) (See recipe in description)

**Step 28.**

leave at room temperature over night (minimum 16 hours)

Eluate the protein

**Step 29.**

Open valve and drip eluate into tubes in 0.5 ml aliquots

**Step 30.**

Measure protein concentration

**Step 31.**

The aliquots with the highest concentrations is mixed and used for further work.

