



Sep 23, 2019

# HiFi Gibson Assembly (Protocol for the NEBuilder® HiFi DNA Assembly Master Mix)

Alba Balletbó<sup>1</sup><sup>1</sup>Wageningen University
1 Works for me
[dx.doi.org/10.17504/protocols.io.7kchksw](https://dx.doi.org/10.17504/protocols.io.7kchksw)

Alba Balletbó ⚡

## ABSTRACT

This is the protocol for DNA Assembly using the NEBuilder® HiFi DNA Assembly Master Mix.

## GUIDELINES

### Optimal Quantities

NEB recommends a total of 0.03–0.2 pmols of DNA fragments when 1 or 2 fragments are being assembled into a vector, and 0.2–0.5 pmols of DNA fragments when 4–6 fragments are being assembled. Efficiency of assembly decreases as the number or length of fragments increases. To calculate the number of pmols of each fragment for optimal assembly, based on fragment length and weight, we recommend the following formula, or using the tool, [NEBcalculator](#).

pmols = (weight in ng) x 1,000 / (base pairs x 650 daltons)

50 ng of 5000 bp dsDNA is about 0.015 pmols

50 ng of 500 bp dsDNA is about 0.15 pmols

The mass of each fragment can be measured using the NanoDrop instrument, absorbance at 260 nm or estimated from agarose gel electrophoresis followed by ethidium bromide staining.

## MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
NEBuilder HiFi DNA Assembly Master Mix - 10 rxns	E2621S	New England Biolabs

- 1 Set up the following reaction on ice (to 20µl total volume):

	2-3 Fragment Assembly	4-6 Fragment Assembly
DNA Ratio	Vector:Insert = 1:2	Vector:Insert = 1:1
Total amount of Fragments	0.03-0.2 pmols	0.2-0.5 pmols
NEBuilder HiFi DNA Assembly Master Mix	10 µL	10 µL
Deionized H2O	10 - X µL	10 - X µL
Total Volume	20 µL	20 µL

- 2 Incubate samples in a thermocycler at 50°C for 15 minutes (when 2 or 3 fragments are being assembled) or 60 minutes (when 4–6 fragments are being assembled).
- 3 Following incubation, store samples on ice or at –20°C for subsequent transformation.

- 4 Transform into chemically competent cells (1-5µL) or in electrocompetent cells (1 µL, diluted 10x or 5 µL, after purification to remove salts).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited