



Oct 14, 2019

Glycerol stocks

[iGEM Dusseldorf](#)¹¹Heinrich-Heine Universität Düsseldorf

1

Works for me

[dx.doi.org/10.17504/protocols.io.76khrw](https://doi.org/10.17504/protocols.io.76khrw)

iGEM Dusseldorf ⚡

ABSTRACT

- Overnight liquid culture (3mL) + Antibiotic.
- For E.coli Low salt LB
- For Pichia YPD
- 80% Glycerol: mix 20mL H₂O with 80mL pure glycerol and then autoclave
- E. coli: 812µL of cells mix with 188 µL of sterile 80% glycerol
- Pichia: Spin 2min at 2000g. Decant 2ml. Resuspend in remaining liquid. 812µL of cells mix with 188µL of sterile 80% glycerol.
- Use cryovials for freezing away.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited