



Aug 12,
2019

NADH Oxidase Activity

Eva Feldman¹

¹University of Michigan - Ann Arbor

1 Works for me dx.doi.org/10.17504/protocols.io.3jcgkiw

Diabetic Complications Consortium
Tech. support email: rmcindoe@augusta.edu

Lili Liang

ABSTRACT

Summary:

Describes assay to quantitate NADH Oxidase activity from tissues.

Diabetic Complications:



Cardiovascular



Nephropathy



Neuropathy



Pediatric Endocrinology



Retinopathy



Uropathy



Wound-Healing

EXTERNAL LINK

<https://www.diacomp.org/shared/document.aspx?id=36&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
50mM Tris	15504-012	Gibco - Thermo Fischer
50mM MES buffer	View	Sigma Aldrich
150µM NADH	N 6879	Sigma Aldrich

MATERIALS TEXT

Reagent Preparation:

Tris-MES buffer (pH 7.0):

NADH

Tris-MES buffer (pH 7.0): Prepare 50mM Tris buffer solution and pH to 7.0 with 50mM MES. 50mM Tris - 302.85mg Tris in 50ml deionized H₂O. 50mM MES - 319.89mg MES in 50ml deionized H₂O.

NADH: Prepare 600µM solution. 2.55mg in 6mL de-ion H₂O, enough for whole plate.

Note:

Fisher Scientific, [RRID:SCR_008452](#)

Sigma-Aldrich, [RRID:SCR_008988](#)

Sample Preparation:

- 1 Turn on Multiskan, set temp to 37°C and set up plate layout.
- 2 Sonicate tissue on ice in 20mM PB pH 7.4 with PMSF inhibitor or thaw prepared samples on ice.
- 3 Remove 25µL for protein analysis.
- 4 Prepare NADH, enough for whole plate.
- 5 Dilute samples 1:5 with de-ionized H₂O.
- 6 Using a clear plate: Add 50µL sample to wells and 50µL diluted sample to wells in duplicate.
- 7 Add 50µL buffer to 3 wells for blanks for positive control.
- 8 Add 100µL Tris-Mes to each sample and blanks.
- 9 For negative control add 200µL Tris-Mes to 3 wells.
- 10 Place plate in Multiskan and add 50µL 600µM NADH to the sample and positive blanks.
Do not add NADH to the 3 negative control wells.
- 11 Press start and read at 340nm for 10 minutes @ 1 minute intervals.
- 12 Save raw data as an Excel file into the NADHx data folder. Use the naming convention NAXXXX.xls, where XXXX is the date in mmdd format.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited