

# RNA extraction from *Synechocystis* sp. PCC 6803 with Trizol reagent

Dennis Dienst

## Abstract

**Citation:** Dennis Dienst RNA extraction from *Synechocystis* sp. PCC 6803 with Trizol reagent. **protocols.io**

dx.doi.org/10.17504/protocols.io.j3scqne

**Published:** 27 Sep 2017

## Protocol

### Step 1.

Label all required tubes and store on ice (or freezer, if not needed immediately)

### Step 2.

Pre-cool all required centrifuges

### Step 3.

Fill sterile 50 mL tube (Falcon) w/ ice and store in ice bath

### Step 4.

Pour 20-25 mL cell culture ( $OD_{750} < 1.0$ ) to ice-filled tube (up to 45 mL mark) Note: avoid long transport of cell culture before harvest

### Step 5.

spin down at 4000 - 5000 g and 4 ° C for 5 min

### Step 6.

discard supernatant (w/ ice) into big beaker Note: depending on strain/mutant some cells will get lost at this step

### Step 7.

resuspend cell pellet in residual water ( 1 mL)

### Step 8.

transfer suspension into 2 mL (safe lock!) tubes (work on ice!)

### Step 9.

spin down at 13.000 g and 4° C for 15 sec

### Step 10.

discard supernatant w/ pipet Note: try to remove supernatant 'as quantitatively as possible'

### Step 11.

resuspend pellet in 1 mL Trizol reagent

### Step 12.

store at -20 °C or (better) -80 °C

### Step 13.

incubate frozen samples at 65° C for 15 min under constant agitation Note: if no shaking thermoblock is available, vortex once per minute

**Step 14.**

add 200 µL (ice cold) chloroform-isoamylalcohol (24:1) per 1 mL Trizol and vortex for 30 sec

**Step 15.**

spin down at 11.000 g and 4° C for 10 min

**Step 16.**

transfer upper, aqueous phase to fresh 1.5 mL tube

**Step 17.**

add 1 Vol. (ice cold) phenol-chloroform-isoamylalcohol (25:24:1) Note: for RNA preparation phenol solution/mixtures should not be Tris-buffered

**Step 18.**

transfer upper, aqueous phase to fresh 1.5 mL tube

**Step 19.**

add 1 Vol. isopropanol (2-propanol), 10 µL 3 M Na-Acetate (pH 5.2) and 1 µL glycogen (RNA grade, Thermo)

**Step 20.**

incubate o/n at -20 for precipitation

**Step 21.**

spin down at 13.000 g and 4° C for 30 min

**Step 22.**

remove supernatant and wash pellet w/ 70% EtOH (ice cold)

**Step 23.**

spin down at 13.000 g and 4° C for 10 min

**Step 24.**

repeat steps 24 and 25

**Step 25.**

discard supernatant and air-dry RNA pellet for 10 min

**Step 26.**

resuspend RNA Pellet in 30 µL ultra-pure water and store at -20 or (better) -80°C