

Cross-linking of IgG to Protein A or G Beads (S1425/S1430)

New England Biolabs

Abstract

This protocol consists of an IgG purification step followed by covalent cross-linking of the IgG to the Protein A/G solid support. For IgG that has been previously purified, proceed directly to the cross-linking protocol.

Citation: New England Biolabs Cross-linking of IgG to Protein A or G Beads (S1425/S1430). [protocols.io](https://doi.org/10.17504/protocols.io.crdv25)

[dx.doi.org/10.17504/protocols.io.crdv25](https://doi.org/10.17504/protocols.io.crdv25)

Published: 10 Feb 2015

Guidelines

Overview

Materials Needed:

Protein A ([NEB #S1425S](#)) or Protein G ([NEB #S1430S](#)) Magnetic Beads

Elution Buffer: 0.1 M glycine-HCl (pH 2.5)

Binding Buffer: 0.1 M NaPhosphate Buffer (pH 8.0)

Dimethyl pimelidate dihydrochloride (Sigma, D-8388) dissolved at 25 mM in Cross-linking Buffer.

Cross-linking Buffer: 0.2 M triethanolamine (pH 8.2)

Blocking Buffer: 0.1 M ethanolamine (pH 8.2)

Immunoglobulin in Binding Buffer

This protocol consists of an IgG purification step followed by covalent cross-linking of the IgG to the Protein A/G solid support. For IgG that has been previously purified, proceed directly to the cross-linking protocol.

Materials

 Protein A Magnetic Beads - 1 ml [S1425S](#) by [New England Biolabs](#)

 Protein G Magnetic Beads - 1 ml [S1430S](#) by [New England Biolabs](#)

Protocol

IgG Purification

Step 1.

Vortex and thoroughly resuspend Protein A Magnetic Beads

IgG Purification

Step 2.

Aliquot **100 µl** of bead suspension to a sterile microcentrifuge tube

 **AMOUNT**

100 µl Additional info:

IgG Purification

Step 3.

(wash #1) Add **500 µl** 0.1 M NaPhosphate Buffer (pH 8.0)

 **AMOUNT**

500 µl Additional info:

IgG Purification

Step 4.

(wash #1) Vortex to resuspend

 **DURATION**

00:30:00

IgG Purification

Step 5.

(wash #1) Apply magnet for 30 seconds, to pull beads to the side of the tube

 **DURATION**

00:00:30

IgG Purification

Step 6.

(wash #1) Remove supernatant

IgG Purification

Step 7.

(wash #2) Add **500 µl** 0.1 M NaPhosphate Buffer (pH 8.0)

 **AMOUNT**

500 µl Additional info:

IgG Purification

Step 8.

(wash #2) Vortex to resuspend

IgG Purification

Step 9.

(wash #2) Apply magnet for 30 seconds, to pull beads to the side of the tube

 **DURATION**

00:00:30

IgG Purification

Step 10.

(wash #2) Remove supernatant

IgG Purification

Step 11.

Add to the beads **80 µl** of 0.1 M NaPhosphate Buffer (pH 8.0)

 **AMOUNT**

80 µl Additional info:

IgG Purification

Step 12.

Add 15-25 µl of serum **OR** 20 µg purified IgG in a maximum volume of 30 µl

IgG Purification

Step 13.

Mix thoroughly and incubate at 4°C with agitation for 30 minutes

 DURATION

00:30:00

IgG Purification

Step 14.

Apply magnet and remove supernatant.

IgG Purification

Step 15.

(wash #1) Add 500 µl 0.1 M NaPhosphate Buffer (pH 8.0)

IgG Purification

Step 16.

(wash #1) Vortex to resuspend

IgG Purification

Step 17.

(wash #1) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Purification

Step 18.

(wash #1) Remove supernatant

IgG Purification

Step 19.

(wash #2) Add 500 µl 0.1 M NaPhosphate Buffer (pH 8.0)

IgG Purification

Step 20.

(wash #2) Vortex to resuspend

IgG Purification

Step 21.

(wash #2) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Purification

Step 22.

(wash #2) Remove supernatant

IgG Purification

Step 23.

(wash #3) Add 500 µl 0.1 M NaPhosphate Buffer (pH 8.0)

IgG Purification

Step 24.

(wash #3) Vortex to resuspend

IgG Purification

Step 25.

(wash #3) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Purification

Step 26.

(wash #3) Remove supernatant

🔗 NOTES

New England Biolabs 03 Feb 2015

At this point the purified IgG can be eluted from the beads or used directly for immunoprecipitation of target proteins. The purified IgG can also be cross-linked to the Protein A beads (see cross-linking protocol) to create a reusable immunoprecipitation bead which prevents the co-elution of antibody with target protein.

IgG Cross-linking

Step 27.

Add **1 ml** of Cross-linking Buffer (0.2 M triethanolamine, [pH 8.2]) to the Protein A/G immobilized antibody (wash **a**)

📄 AMOUNT

1 ml Additional info:

🔗 NOTES

New England Biolabs 03 Feb 2015

At this point the purified IgG can be eluted from the beads or used directly for immunoprecipitation of target proteins. The purified IgG can also be cross-linked to the Protein A beads (see cross-linking protocol) to create a reusable immunoprecipitation bead which prevents the co-elution of antibody with target protein.

IgG Cross-linking

Step 28.

Vortex to resuspend (wash **a**)

IgG Cross-linking

Step 29.

Apply magnet for 30 seconds, to pull beads to the side of the tube (wash **a**)

🕒 DURATION

00:00:30

IgG Cross-linking

Step 30.

Remove supernatant (wash **a**)

IgG Cross-linking

Step 31.

Add 1 ml of Cross-linking Buffer (0.2 M triethanolamine, [pH 8.2]) to the Protein A/G immobilized antibody (wash **b**)

IgG Cross-linking

Step 32.

Vortex to resuspend (wash **b**)

IgG Cross-linking

Step 33.

Apply magnet for 30 seconds, to pull beads to the side of the tube (wash **b**)

🕒 DURATION

00:00:30

IgG Cross-linking

Step 34.

Remove supernatant (wash **b**)

IgG Cross-linking

Step 35.

Resuspend in 1 ml Cross-linking Buffer containing 25 mM DMP (6.5 mg DMP/ml of buffer)

IgG Cross-linking

Step 36.

Mix thoroughly and incubate at room temperature for 45 minutes with agitation

 DURATION

00:45:00

IgG Cross-linking

Step 37.

Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 38.

Remove supernatant

IgG Cross-linking

Step 39.

Add 1 ml Blocking Buffer (0.1 M ethanolamine, [pH 8.2])

IgG Cross-linking

Step 40.

Vortex to resuspend

IgG Cross-linking

Step 41.

Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 42.

Remove supernatant

IgG Cross-linking

Step 43.

Add 1 ml of Blocking Buffer

 AMOUNT

1 ml Additional info:

IgG Cross-linking

Step 44.

Vortex to resuspend

IgG Cross-linking

Step 45.

Incubate for 1 hour at room temperature with agitation

 DURATION

01:00:00

IgG Cross-linking

Step 46.

Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 47.

Remove supernatant

IgG Cross-linking

Step 48.

(wash #1) Add 1 ml of PBS

 AMOUNT

1 ml Additional info:

IgG Cross-linking

Step 49.

(wash #1) Vortex to resuspend

 AMOUNT

1 ml Additional info:

IgG Cross-linking

Step 50.

(wash #1) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 51.

(wash #1) Remove supernatant

IgG Cross-linking

Step 52.

(wash #2) Add 1 ml of PBS

IgG Cross-linking

Step 53.

(wash #2) Vortex to resuspend

IgG Cross-linking

Step 54.

(wash #2) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 55.

(wash #2) Remove supernatant

IgG Cross-linking

Step 56.

(wash #3) Add 1 ml of PBS

IgG Cross-linking

Step 57.

(wash #3) Vortex to resuspend

IgG Cross-linking

Step 58.

(wash #3) Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

IgG Cross-linking

Step 59.

(wash #3) Remove supernatant

IgG Cross-linking

Step 60.

Add 1 ml Elution Buffer (0.1 M glycine-HCl [pH 2.5])

 AMOUNT

1 ml Additional info:

IgG Cross-linking

Step 61.

Vortex to resuspend

IgG Cross-linking

Step 62.

Apply magnet for 30 seconds, to pull beads to the side of the tube

 DURATION

00:00:30

 NOTES

New England Biolabs 02 Feb 2015

This elutes bound antibody that is not cross-linked with DMP.

IgG Cross-linking

Step 63.

Remove supernatant

IgG Cross-linking

Step 64.

Resuspend and store beads in 100 µl PBS, 0.1% Tween 20, 0.02% sodium azide