



Apr 15, 2019

Working

Centrifuge the whole blood to separate red blood cells

[George Lykotrafitis](#)¹¹University of Connecticut[dx.doi.org/10.17504/protocols.io.zzwf77e](https://doi.org/10.17504/protocols.io.zzwf77e) [George Lykotrafitis](#) 

ABSTRACT

Centrifuge the whole blood in order to obtain the red blood cells.

MATERIALS TEXT

Human whole blood, centrifuge, centrifuge tubes.

- 1 Centrifuge the whole blood using the parameters: 500 ×g, 4 °C, and 5 min.
- 2 Remove the plasma and buffy coat.
- 3 Add the same volume of Alsever's solution as the remaining blood components.
- 4 Centrifuge again using the parameters: 500 ×g, 4 °C, and 5 min.
- 5 Remove the top layer.
- 6 Add the same volume of Alsever's solution as the remaining blood components.
- 7 Centrifuge a third time using the parameters: 500 ×g, 4 °C, and 5 min.
- 8 Remove the top layer.
- 9 Add the same volume of Alsever's solution as the remaining blood component - red blood cells, and store at 4 °C for up to 7 days.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited