

Electroporation of *Oxyrrhis marina*

Yutaka Hanawa¹, Yoshihisa Hirakawa¹, Patrick Keeling²

¹University of Tsukuba, ²University of British Columbia

dx.doi.org/10.17504/protocols.io.vcne2ve

Protist Research to Optimize Tools in Genetics (PROT-G)




Yoshihisa Hirakawa

University of Tsukuba



ABSTRACT



PROTOCOL STATUS

In development

We are still developing and optimizing this protocol

Culture condition of *O. marina*

- 1 Transfer *Oxyrrhis marina* (NIES-494) cells to 20 mL of fresh IMK medium (Nihon Pharmaceutical Co., Ltd.) in a plastic flask (IWAKI 75 cm²) at the concentration of 200 cells/mL, and add *Pyramimonas parkeae* (NIES-254) as feed at the concentration of 1×10^4 cells/mL.

Grow cells at 22°C with a light/dark cycle of 14h/8h for 14 days. After two weeks, the cell density of *O. marina* will reach approximately 1×10^4 cells/mL.

Concentration of *O. marina* cells

- 2 Collect *O. marina* and *P. parkeae* cells from 20 mL culture by centrifugation at 1,500 g for 4 min with a swing rotor.

Resuspend cells with 2 mL fresh IMK medium, and then transfer to a 12-well plastic plate. Incubate the plate in a dark condition at 22°C.

After 4 days, the cell density of *O. marina* will be increased from 1×10^5 cells/mL (0 day) to 5×10^5 cells/mL (4 days), in contrast with the drastic decrease of *P. parkeae* cells.

Electroporation with Bio-Rad Gene Pulser Xcell

- 3 Harvest *O. marina* cells at 1,500 g for 4 min with a swing rotor, and then resuspend the cell pellet by 100 μ L Gene Pulser electroporation buffer (Bio-Rad #1652676) at the final concentration of 1×10^6 to 5×10^6 cells/mL.

Add DNA (5-25 μ g) or RNA (5 μ g) to the cell solution, and transfer it to a 0.2 cm cuvette (Bio-Rad). Immediately after electrophoresis, add 1.9 mL fresh IMK medium into the cuvette, and transfer the cells to a 12-well plastic plate.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited