

untitled protocol

Yuan Yao


Abstract

Citation: Yuan Yao untitled protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.qabdsan


Published: 22 May 2018


Materials

 Nitrophenyl b-d-galactopyranoside view by P212121

 Magnesium Sulfate (MgSO₄) Solution - 6.0 ml B1003S by New England Biolabs

✓ Sodium Phosphate monobasic by Contributed by users

 Sodium phosphate dibasic 7558-79-4 by Sigma Aldrich

 500g Sodium Carbonate; Na₂CO₃ (anhydrous) RC-125 by G-Biosciences

✓ KCl by Contributed by users

 1% β-mercaptoethanol by BBI Biotech

Protocol

Exponentially growing cells (1 ml) at 37°C in ABTGcasa medium were collected at OD₄₅₀=0.1, 0.2, 0.3, 0.4 and 0.5

Step 1.

Mix with cold toluene (0.1 ml) and kept on ice immediately

Step 2.

0.2 ml toluene-treated sample was added to 1 ml Z buffer (40 mM NaH₂PO₄, 60 mM Na₂HPO₄, 10 mM KCl, 1 mM MgSO₄ and 50 mM β-mercaptoethanol, pH 7.0) containing 0.66 mg/ml o-nitrophenyl-β-D-galactopyranoside.

Step 3.

The reaction was performed at 30°C until the color changed to yellow

Step 4.

Stop by addition of 0.5 ml 1 M Na₂CO₃

Step 5.

The absorbance at OD₄₂₀ was measured

Step 6.

The β-galactosidase activity was calculated by $1000 \times \text{OD}_{420} / \text{reaction time (min)} \times \text{OD}_{450} \times 0.2 \text{ ml}$

Step 7.