

Basic components detection of Dendrobium plants

Cheng Song

Abstract

Citation: Cheng Song Basic components detection of Dendrobium plants. **protocols.io**

dx.doi.org/10.17504/protocols.io.h4gb8tw

Published: 19 May 2017

Protocol

Extraction

Step 1.

- 1.Plants collected were washed and dried by filter paper.
- 2.Fixed at 105 °C for 20 min and dried to a constant weight at 60 °C.
- 3.One gram of powder with 20 mL of methanol was extracted twice using an ultrasonic cleaner at ambient temperature.
- 4.The extraction solution was centrifuged for 10 min at 5000×g and supernatant was evaporated to dryness.
- 5.Dissolved with 2 mL of dilute sulfuric acid (0.1%) and filtered through a 0.22 µm nylon syringe.

Purification

Step 2.

The extraction was pretreated using Oasis MCX SPE cartridges.The procedure was as follow.

- 1.activating with 2 mL methanol
- 2.conditioning with 2 mL pure water
- 3.loading with 2mL sample solution
- 4.washing with 2 mL formic acid (2%)
- 5.washing with 2 mL methanol
- 6.eluting with 2 mL ammonia methanol (5%)

7.drying with nitrogen, and reconstituting with 2 mL formic acid (0.1%)

HPLC analysis

Step 3.

1.Instruments & materials

Agilent 1260 series HPLC instrument

Xcalibur 2.1 (Thermo Fisher Scientific)

Waters Atlantis T3 column (150 mm × 4.6 mm id, 3 μm)

2.Methods

Mobile phases: acetonitrile (A) with 0.1% formic acid and ultrapure water (B) with 0.1% formic acid.

Gradient elution process: 0-2 min (0% A), 2-12 min (0 to 15% A), 12-22 min (15 to 35% A), 22-32 min (35 to 80% A), and 32-37 min (80 to 0% A).

volume flow: 1.2 mL/min

detection wavelength: 280 nm

column temperature: 25 centigrade

injection volume 10 μl

Mass spectra conditions

Step 4.

1.The ion-source parameters were as follows: sheath gas at 25 arb, auxiliary gas at 3 arb, spray voltage at 4 kV, capillary temperature at 320 °C, tube lens at 120 V, and capillary voltage at 30 V.

2.MS data were collected at $100 \leq m/z \leq 1000$ in positive-ion mode.

3.Two data acquisition methods were used in the experiments. In one method, a high-resolution scan was conducted using the Orbitrap mass analyzer to acquire the MS data at a resolution at 30,000 FWHM and the MS² data at 15,000 FWHM. In the other method, a high-resolution scan was conducted using the Orbitrap to acquire the MS data at a resolution at 30,000 FWHM and the LTQ dynode was used for scanning the MS² and MS³ spectra.

4.Collision-induced dissociation (CID) energy was set at 40%.

Mass spectra interpretation

Step 5.

1. Peaks could be extracted from total ion current chromatograms manually.
2. Or ran XCMS package under R environment to get peak lists.
3. Obtain putative molecular formulae by isotopic patterns verification and the accurate molecular weight.
4. Matching Metlin, Dictionary of Natural Products and ChempSpider databases to identify tentatively.
5. Using Mass Frontier software to speculate the fragmentation pathway for compounds without MSⁿ spectra.